# FUNDACION Dr.JORDI MAS

REVISAT 2008

## ETHICON Wound Closure Manual

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## **Referencias**

1. Taylor EJ et al (ed): *Dorland's Medical Dictionary*, Philadelphia: Saunders, 1988
2. Davis JH: *Clinical Suregery*, CV Mosby, 1987
3. American College of Surgeons, Altemeir WA, et al (eds): *Manual on Control of Infection in Surgical Patients*, Philadelphia: Lippincott, 1976
4. Cruse P, Foorde R: The epidemology of wound infencion, *Surg Clin NA* 60(1):27-40, 1980
5. Gottrup F: *Acta Chir Scand* 149:57, 1987
6. Chu C, Williams DF: Effects of physical configuration and chemical structure of suture materials on bacterial adhesion, *Am J Surg* 147:197-204, 1984
7. Kenady DE: Management of abdominal wounds. *Surg Clin NA* 64(4):803-806, 1984

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SUTURE	TYPE	COLOR OF MATERIAL	RAW MATERIAL
Surgical Gut Suture	Plain	Yellowish-tan Blue Dyed	Collagen derived from serosa of beef intestine or submucosa of sheep intestine.
Surgical Gut Suture	Chromic	Brown Blue Dyed	Collagen derived from serosa of beef intestine or submucosa of sheep intestine.
Coated VICRYL RAPIDE (polyglactin 910) Suture	Braided	Undyed (Natural)	Copolymer of glycolide and lactide coated with polyglactin 370 and calcium stearate.
Coated VICRYL (polyglactin 910) Suture	Braided Monofilament	Violet Undyed (Natural)	Copolymer of glycolide and lactide coated with polyglactin 370 and calcium stearate.
MONOCRYL (poliglecaprone 25) Suture	Monofilament	Undyed (Natural) Violet	Copolymer of glycolide and epsilon-caprolactone.
PDS II (polydioxanone) Suture	Monofilament	Violet Blue Clear	Polyester polymer.
PANACRYL Suture	Braided	Undyed (Natural)	Copolymer of lactide and glycolide coated with a polymer of caprolactone and glycolide.
PERMA-HAND Silk Suture	Braided	Black White	Organic protein called fibroin.
Surgical Stainless Steel Suture	Monofilament Multifilament	Silver metallic	316L stainless steel.
ETHILON Nylon Suture	Monofilament	Black Green Undyed (Clear)	Long chain aliphatic polymers Nylon 6 or Nylon 6,6.
NUROLON Nylon Suture	Braided	Black Green Undyed (Clear)	Long-chain aliphatic polymers Nylon 6 or Nylon 6,6.
MERSILENE Polyester Fiber Suture	Braided Monofilament	Green Undyed (White)	Poly (ethylene terephthalate).
ETHIBOND EXCEL Polyester Fiber Suture	Braided	Green Undyed (White)	Poly (ethylene terephthalate) coated with polybutylate.
PROLENE Polypropylene	Monofilament	Clear Blue	Isotactic crystalline stereoisomer of polypropylene.
PRONOVA Poly(vinylidene fluoride)	Monofilament	Clear Blue	Polymer blend of poly(vinylidene fluoride) and

propylene-VDF) Suture			poly (vinylidene fluoride- cohexafluoropropylene)
-----------------------	--	--	--

Suturing Options:  
Materials, Characteristics and Applications Table  
Part 2

SUTURE	TENSILE STRENGTH RETENTION <i>IN VIVO</i>	ABSORPTION RATE	TISSUE REACTION	CONTRAINDICATIONS
Surgical Gut Suture	Individual patient characteristics can affect rate of tensile strength loss.	Absorbed by proteolytic enzymatic digestive process.	Moderate reaction	Being absorbable, should not be used where extended approximation of tissues under stress is required. Should not be used in patients with known sensitivities or allergies to collagen or chromium.
Surgical Gut Suture	Individual patient characteristics can affect rate of tensile strength loss.	Absorbed by proteolytic enzymatic digestive process.	Moderate reaction	Being absorbable, should not be used where extended approximation of tissues under stress is required. Should not be used in patients with known sensitivities or allergies to collagen or chromium.
Coated VICRYL RAPIDE (polyglactin 910) Suture	Approximately 50% remains at 5 days. All tensile strength is lost by approximately 10-14 days.	Essentially complete by 42 days. Absorbed by hydrolysis.	Minimal to moderate acute inflammatory reaction.	Should not be used when extended approximation of tissue under stress is required or where wound support beyond 7 days is required.
Coated VICRYL (Polyglactin 910) Suture	Approximately 75% remains at 2 weeks. Approximately 50% remains at 3 weeks. (6-0 and larger)	Essentially complete between 56-70 days. Absorbed by hydrolysis.	Minimal acute inflammatory reaction	Being absorbable, should not be used where extended approximation of tissue is required.
MONOCRYL (poliglecaprone 25) Suture	Approximately 50-60% remains at 1 week (undyed); 60-70% at 1 week (dyed). Approximately 20-30% remains at 2 weeks (undyed); 30-40% at 2 weeks (dyed). Lost by 28 days.	Complete at 91-119 days. Absorbed by hydrolysis.	Minimal acute inflammatory reaction	Being absorbable, should not be used where extended approximation of tissue under stress is required. Undyed not indicated for use in fascia.
PDS II (polydioxanone) Suture	Approximately 70% remains at 2 weeks. Approximately 50% remains at 4 weeks. Approximately 25%	Minimal until about 90th day. Essentially complete within 6 months	Slight reaction	Being absorbable, should not be used where prolonged approximation of tissues under stress is required. Should not be

	remains at 6 weeks.	Absorbed by slow hydrolysis.		used with prosthetic devices, such as heart valves or synthetic grafts.
PANACRYL Suture	Approximately 80% remains at 3 months.Approximately 60% remains at 6 months.	Essentially complete in 1.5 to 2.5 years.	Minimal inflammatory reaction	Being absorbable, should not be used where extended approximation of tissue beyond 6 months is required.
PERMA-HAND Silk Suture	Progressive degradation of fiber may result in gradual loss of tensile strength over time.	Gradual encapsulation by fibrous connective tissue.	Acute inflammatory reaction	Should not be used in patients with known sensitivities or allergies to silk.
Surgical Stainless Steel Suture	Indefinite.	Nonabsorbable.	Minimal acute inflammatory reaction	Should not be used in patients with known sensitivities or allergies to 316L stainless steel, or constituent metals such as chromium and nickel.
ETHILON Nylon suture	Progressive hydrolysis may result in gradual loss of tensile strength over time.	Gradual encapsulation by fibrous connective tissue.	Minimal acute inflammatory reaction	Should not be used where permanent retention of tensile strength is required.
NUROLON Nylon Suture	Progressive hydrolysis may result in gradual loss of tensile strength over time.	Gradual encapsulation by fibrous connective tissue.	Minimal acute inflammatory reaction	Should not be used where permanent retention of tensile strength is required.
MERSILENE Polyester Fiber Suture	No significant change known to occur in vivo.	Gradual encapsulation by fibrous connective tissue.	Minimal acute inflammatory reaction	None known.
ETHIBOND EXCEL Polyester Fiber Suture	No significant change known to occur in vivo.	Gradual encapsulation by fibrous connective tissue.	Minimal acute inflammatory reaction	None known.
PROLENE Polypropylene	Not subject to degradation or weakening by action of tissue enzymes.	Gradual encapsulation by fibrous connective tissue.	Minimal acute inflammatory reaction	None known.

PRONOVA Poly (hexafluoropropylene- VDF) Suture	Not subject to degradation or weakening by action of tissue enzymes.	Gradual encapsulation by fibrous connective tissue.	Minimal to mild inflammatory reaction	None known.
---	---	---	--	-------------

Suturing Options:  
Materials, Characteristics and Applications Table  
Part 3

SUTURE	FREQUENT USES	HOW SUPPLIED	COLOR CODE
Surgical Gut Suture	General soft tissue approximation and/or ligation, including use in ophthalmic procedures. Not for use in cardiovascular and neurological tissues.	7-0 thru 3 with and without needles, and on LIGAPAK dispensing reels 0 thru 1 with CONTROL RELEASE needles	Yellow
Surgical Gut Suture	General soft tissue approximation and/or ligation, including use in ophthalmic procedures. Not for use in cardiovascular and neurological tissues.	7-0 thru 3 with and without needles, and on LIGAPAK dispensing reels 0 thru 1 with CONTROL RELEASE needles	Beige
Coated VICRYL RAPIDE (polyglactin 910) Suture	Superficial soft tissue approximation of skin and mucosa only. Not for use in ligation, ophthalmic, cardiovascular, or neurological procedures.	5-0 thru 1 with needles	Violet and Red
Coated VICRYL (Polyglactin 910) Suture	General soft tissue approximation and/or ligation, including use in ophthalmic procedures. Not for use in cardiovascular and neurological tissues.	8-0 thru 3 with and without needles, and on LIGAPAK dispensing reels 4-0 thru 2 with CONTROL RELEASE needles 8-0 with attached beads for ophthalmic use	Violet
MONOCRYL (poliglecaprone 25) Suture	General soft tissue approximation and/or ligation. Not for use in cardiovascular and neurological tissues, microsurgery, or ophthalmic surgery.	6-0 thru 2 with and without needles 3-0 thru 1 with CONTROL RELEASE needles	Coral
PDS II (polydioxanone) Suture	All types of soft tissue approximation, including pediatric cardiovascular and ophthalmic procedures. Not for use in adult cardiovascular tissue, microsurgery and neural tissue.	9-0 thru 2 with needles (violet) 4-0 thru 1 with CONTROL RELEASE needles (violet) 9-0 thru 7-0 with needles (blue) 7-0 thru 1 with needles (clear)	Silver
PANACRYL Suture	General soft tissue approximation and/or ligation. Not for use in ophthalmic, cardiovascular or neurological	3-0 thru 2 with CONTROL RELEASE needles or without needles	Plum

	tissues.		
PERMA-HAND Silk Suture	General soft tissue approximation and/or ligation, including cardiovascular, ophthalmic and nuerological procedures.	9-0 thru 5 with and without needles, and on LIGAPAK dispensing reels 4-0 thru 1 with CONTROL RELEASE needles	Light Blue
Surgical Stainless Steel Suture	Abdominal wound closure, hernia repair, sternal closure, and orthopaedic procedures including cerclage and tendon repair.	10-0 thru 7 with and without needles	Yellow-Ochre
ETHILON Nylon suture	General soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.	11-0 thru 2 with and without needles	Mint Green
NUROLON Nylon Suture	General soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.	6-0 thru 1 with and without needles 4-0 thru 1 with CONTROL RELEASE needles	Mint Green
MERSILENE Polyester Fiber Suture	General soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.	6-0 thru 5 with or without needles 10-0 and 11-0 for ophthalmic (green monofilament) 0 with CONTROL RELEASE needles	Turquoise
ETHIBOND EXCEL Polyester Fiber Suture	General soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.	7-0 thru 5 with or without needles 4-0 thru 1 with CONTROL RELEASE needles-various sizes attached to TFE polymer pledgets	Orange
PROLENE Polypropylene Suture	General soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.	7-0 thru 2(clear) with and without needles 10-0 thru 8-0 and 6-0 thru 2(blue) with and without needles 0 thru 2 with CONTROL RELEASE needles-various sizes attached to TFE polymer pledgets	Deep Blue
PRONOVA Poly (hexafluoropropylene-VDF)Suture	General soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.	6-0 thru 2 (clear) with needles 10-0 thru 8-0 and 6-0 thru 2 (blue) with needles 7-0 are U.S.P. except for diameter 0 thru 2 with	Cranberry

		CONTROL RELEASE needles	
--	--	----------------------------	--

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\*Marca registrada

## Endoscopic Knot Tying Techniques

During an endoscopic procedure, a square knot or surgeon's knot may be tied either outside the abdomen and pushed down into the body through a trocar (extracorporeal) or directly within the abdominal cavity (intracorporeal).

In *extracorporeal knot tying*, the suture appropriately penetrates the tissue, and both needle and suture are removed from the body cavity, bringing both suture ends outside of the trocar. Then a series of half-hitches are tied, each one being pushed down into the cavity and tightened with an endoscopic knot pusher.

*Intracorporeal knot tying* is performed totally within the abdominal cavity. After the suture has penetrated the tissue, the needle is cut from the suture and removed. Several loops are made with the suture around the needleholder, and the end of the suture is pulled through the loops. This technique is then repeated to form a surgeon's knot, which is tightened by the knot pusher.

In both extracorporeal and intracorporeal knot tying, the following principles of suture manipulation on tissue should be observed.

1. Handle tissue as gently as possible to avoid tissue trauma.
2. Grasp as little tissue as possible.
3. Use the smallest suture possible for the task.
4. Exercise care in approximating the knot so that the tissue being approximated is not strangulated.
5. Suture must be handled with care to avoid damage.

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## Referencias

- VanWinkle Jr W, Hastings JC: Considerations in the choice of suture material for various tissues, *Surg Gyn & Obs* 135:113-126, July 1972
- Dineen P: The effect of suture material in the development of vascular infection. *Vasc Surg* 11(1):29-32, January/February 1977
- Blaydes JE, Berry J: Comparative evaluation of coated and uncoated polyglactin 910 in cataract and muscle surgery, *Ophth Surg* 11(11):790-793, November 1980
- Murray DH, Blaisdell FW: Use of synthetic absorbable sutures for abdominal and chest wound closure, *Arch Surg* 113:477-480, April 1978
- Hunter P: 10-0 Vicryl suture offers alternative to no stitch, *Ocular Surg News* 9(16): August 15, 1991
- Craig PH, et al: A biologic comparison of polyglactin 910 and polyglycolic acid synthetic absorbable sutures, *Surg Gyn & Obs* 14(1):1-10, July 1975
- Amshel AL: The use of Vicryl (polyglactin 910) sutures in colonic and rectal surgery, *Dis Col & Rect* 20(7):635-638, October 1977
- Chu CC, Williams DF: Effects of physical configuration and chemical structure of suture material on bacterial adhesion, *Am J Surg* 147:197-204, February 1984
- Pfeffer WH: How to choose sutures for microsurgery, *Contemporary OB/GYN*:1-5, November 1983
- Chusak RB, Dibbell DG: Clinical experience with polydioxanone monofilament absorbable sutures in plastic surgery, *Plas & Recon Surg* 72(2):217-220, August 1983
- Houdart R et al: Polydioxanone in digestive surgery, *Am J Surg* 152:268-271, September 1986
- Durley P, Bucknall TE: Assessment of sutures for use in colonic surgery: An experimental study, *J Royal Society of Med* 77:472-477, June 1984
- Lo Piccolo M (ed): Mersilene suture results found comparable to nylon, *Ophthalmology Times*:18, August 15, 1988
- Kronenthal RL: Intraocular degradation of nonabsorbable sutures, *Amer Intra-Ocular Implant Soc J* 3:222-228, July/October 1977
- Nealon TF: *Fundamental Skills in Surgery*, 3d ed, Philadelphia: Saunders, 1979, p. 47
- Herrman JB: Tensile strength and knot security of surgical suture materials, *Am Surg* 37:211, April 1971
- VanWinkle Jr W: Epithelization and wound contraction, In *Symposium on Basic Science in Plastic Surgery* (Krizek TJ, Hoopes JE, eds). Mosby, St. Louis 1976, p. 87
- Quinn J, et al: A randomized trial comparing octylcyanoacrylate tissue adhesive and sutures in the management of lacerations, *J Am Med Assn* 277(19):1527-1530, May 21, 1997
- Toriumi DM, et al: Use of Octyl-2-Cyanoacrylate for Skin Closure in Facial Plastic Surgery, *Plas & Recon Surg* 102(6):2209-2219, November 1998
- Vallors B, et al: Absorbable or Nonabsorbable Suture Materials for Closure of the Dura Mater?, *Neurosurgery* 9(4):407-413, 1981
- Tera H, Aberg C: Tissue holding power to a single suture in different parts of the alimentary tract, *Acta Chir Scand* 142:343-348, 1976
- McCollum CN, Kester RC: A sliding suture technique for inaccessible arterial anastomoses, *Surg Gynecol Obstet* 153(6):907, 1981
- VanWinkle Jr W: The healing of skin and visceral wounds and the effect of suture materials, In *Symposium on Basic Science in Plastic Surgery* (Krizek TJ, Hoopes JE, eds). Mosby, St. Louis 1976, p. 135

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# ETHICON needle codes and their meaning

## Part 1

CODE	MEANING
<b>BB</b>	Blue Baby
<b>BIF</b>	Blunt Intraocular Fixation
<b>BN</b>	Bunnell
<b>BP</b>	Blunt Point
<b>BV</b>	Blood Vessel
<b>BVH</b>	Blood Vessel Half
<b>C</b>	Cardiovascular
<b>CC</b>	Calcified Coronary
<b>CCS</b>	Conventional Cutting Sternotomy
<b>CE</b>	Cutting Edge
<b>CFS</b>	Conventional for Skin
<b>CIF</b>	Cutting Intraocular Fixation
<b>CP</b>	Cutting Point
<b>CPS</b>	Conventional Plastic Surgery
<b>CPX</b>	Cutting Point Extra Large
<b>CS</b>	Corneal-Scleral
<b>CSB</b>	Corneal-Scleral Bi-Curve
<b>CSC</b>	Corneal-Scleral Compound Curve
<b>CT</b>	Circle Taper
<b>CTB</b>	Circle Taper Blunt
<b>CTX</b>	Circle Taper Extra Large
<b>CTXB</b>	Circle Taper Extra Large Blunt
<b>CV</b>	Cardiovascular
<b>DC</b>	Dura Closure
<b>DP</b>	Double Point
<b>EN</b>	Endoscopic Needle

CODE	MEANING
<b>FSLX</b>	For Skin Large Extra
<b>G</b>	Greishaber
<b>GS</b>	Greishaber Spatula
<b>J</b>	Conjunctive
<b>KS</b>	Keith Straight
<b>LH</b>	Large Half (circle)
<b>LR</b>	Large Retention
<b>LS</b>	Large Sternotomy
<b>M</b>	Muscle
<b>MF</b>	Modified Ferguson
<b>MH</b>	Medium Half (circle)
<b>MO</b>	Mayo
<b>MOB</b>	Mayo Blunt
<b>OPS</b>	Ocular Plastic Surgery
<b>OS</b>	Orthopaedic Surgery
<b>P</b>	Plastic
<b>PC</b>	Precision Cosmetic
<b>PS</b>	Plastic Surgery
<b>RB</b>	Renal (artery) Bypass
<b>RD</b>	Retinal Detachment
<b>RH</b>	Round Half (circle)
<b>RV</b>	Retinal-Vitreous
<b>S</b>	Spatula
<b>SC</b>	Straight Cutting
<b>SFS</b>	Spatulated for Skin
<b>SH</b>	Small Half (circle)

<b>SIF</b>	Ski Intraocular Fixation
<b>SKS</b>	Sternotomy Keith Straight
<b>SM</b>	Spatulated Muscle
<b>ST</b>	Straight Taper

ETHICON needle codes and their meaning  
Part 2

<b>CODE</b>	<b>MEANING</b>
<b>STB</b>	Straight Blunt
<b>STC</b>	Straight Cutting
<b>STP</b>	Straight Taper Point
<b>TE</b>	Three-Eighths
<b>TF</b>	Tetralogy of Fallot
<b>TG</b>	Transverse Ground
<b>TGW</b>	Tranverse Ground Wire
<b>TN</b>	Trocar Needle
<b>TP</b>	Taper Pericostal/Point
<b>TPB</b>	Taper Pericostal/Point Blunt
<b>TS</b>	Tendon Straight
<b>TQ</b>	Twisty Q
<b>UCL</b>	5/8 Circle Colateral Ligament
<b>UR</b>	Urology
<b>URB</b>	Urology Blunt
<b>V</b>	TAPERCUT Surgical Needle
<b>VAS</b>	Vas Deferens
<b>X or P</b>	Exodontal (dental)
<b>XLH</b>	Extra Large Half (circle)
<b>XXLH</b>	Extra Extra Large Half (circle)

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## **References**

- Edlich RF et al: Scientific basis for selecting surgical needles and needle holders for wound closure, *Clin Plastic Surg* 17(3):583-601, July 1990
- Thacker JG et al: Surgical needle sharpness, *Am J Surg* 157:334-339, March 1989
- Trier WC: Considerations in the choice of surgical needles, *Surg Gynecol Obstet* 149:84, 1979
- Towler MA et al: Laser-drilled cardiovascular surgical needles, *J Appl Biomaterials* 2:183-186, 1991
- Troutman RC, McGregor W: New compound curved needle TG6-C Plus, *Ophthalmology, Instrument and Book Suppl* 92(8):84-86, 1985
- Abidin MR et al: Metallurgically bonded needle holder jaws, *Am Surg* 56(10):643-647, October 1990

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## Primary Packets Page 1

Individual sutures and multiple suture strands are supplied sterile within a primary packet. The exterior surfaces of the overwrap are not sterile. ETHICON primary packaging is designed to permit fast and easy opening in one peelable motion. The single layer overwrap of primary packaging is made of either foil *or* coated TYVEK on one side heat-sealed to polyethylene film on the other. Absorbable sutures are always encased in foil to provide a safe and durable moisture barrier and to withstand sterilization in the manufacturing process. Most nonabsorbable sutures are encased in coated TYVEK overwraps. In a continuous effort to be more environmentally conscious, ETHICON has chosen materials in the manufacture of primary packets which generate minimal negative impact to the environment upon incineration or disposal. Furthermore, wherever possible, the number of primary packaging layers has been reduced by as much as 50 percent, thus reducing the volume of environmental water per OR procedure. Each primary packet provides critical product information and the same color coding as its dispenser box. The packet also identifies the product code number, material, size, needle type, and the number of needles per packet to simplify needle counts.

		Table 1
Absorbable Sutures	Layers	Most Common ETHICON Suture Packaging
Surgical Gut Suture (plain and Chromic)	TYVEK overwrap, foil primary package containing one-step RELAY suture delivery system tray	
Coated VICRYL RAPIDE (polyglactin 910) Suture	TYVEK overwrap, foil primary package, paper folder	
Coated VICRYL (polyglactin 910) Suture	Peelable foil overwrap, one-step RELAY tray	
MONOCRYL* (poliglecaprone 25) Suture	Peelable foil overwrap, one-step RELAY tray	
PDS II (polydioxanone) Suture	Peelable foil overwrap, one-step RELAY tray	
Panacryl Suture	Peelable foil overwrap, one-step RELAY tray	
Non-Absorbable Sutures	Layers	
PERMA-HAND Silk Suture	TYVEK overwrap, one-step RELAY tray	
Stainless Steel Suture	TYVEK overwrap, paper folder	
ETHILON Nylon Suture	Peelable foil overwrap, one-step RELAY tray	
NUROLON Nylon Suture	TYVEK overwrap, one-step RELAY tray	
MERSILENE	TYVEK overwrap, one-step RELAY tray	
ETHIBOND EXCEL Polyester Fiber Suture	TYVEK overwrap, one-step RELAY tray	

PROLENE* Polypropylene Suture	TYVEK overwrap, one-step RELAY tray
PRONOVA Poly (hexafluoropropylene-VDF) Suture	TYVEK overwrap, one-step RELAY tray


Primary packets of suture material may contain sutures in one of five styles:

1. Standard lengths of nonneedled material: 54 inches (135cm) of absorbable or 60 inches (150cm) of nonabsorbable suture, which may be cut in half, third, or quarter lengths for ligating or threading.
2. SUTUPAK pre-cut sterile suture is nonneedled material for ligating or threading. These lengths may be supplied in a multistrand labyrinth packet or in a folder packet, both of which are designed to deliver one strand at a time. SUTUPAK sutures may be removed from the packet and placed in the suture book.
3. One single strand of material with single- or double-armed swaged needle(s). Needles for one-step RELAY suture packets, microsurgery, and some ophthalmic needles are secured in a "needle park." The needle park is designed to provide a standard location for, and easy access to, the needle. All other needles are protected within an inner folder or other specific channel within a paper folder.

Most single strand needled sutures are sealed in convenient one-step RELAY delivery packages. One-step RELAY packages allow the needle to be armed in the needleholder from any angle without touching the needle. This increases the safety of handling needles intraoperatively. If it is preferred to locate the needle by hand, this can be accomplished with the one-step RELAY package by pushing up the flap behind the needle park, thereby elevating the needle so it can be grasped by hand.

**FIG.2**

**Method for  
Preparing One-  
Step Relay  
Package  
Sutures**

Arm the needle directly from the one-step RELAY tray and deliver the single suture to the surgeon.

## Primary Packets Page 2

Multiple suture strands, either swaged to a single needle or double-armed. This type is appropriate for procedures requiring numerous interrupted sutures of the same type. It saves valuable operative time by enabling the surgeon to use one suture while the next is being armed--without delay of opening packets or threading needles.

Multistrand packets are labeled with the symbol MS/ that denotes multiple strands/number of strands of surgical needles per packet. Multistrand packets may contain 3 to 10 swaged sutures. The inner folder for these products is white.

All packets containing CONTROL RELEASE needle sutures have multiple strands (8, 5, 4, 3, or 1) and are designated CR/8, CR/5, CR/4, CR/3, or CR/1. CONTROL RELEASE sutures may be available in foil or Tyvek overwrap packets for single strand delivery. The single strand delivery *folder* is used for some Coated VICRYL (polyglactin 910) sutures, MONOCRYL (poliglecaprone 25) sutures, PDS II (polydioxanone) sutures, ETHIBOND EXCEL polyester sutures, NUROLON nylon sutures, MERSILENE polyester sutures, and PERMA-HAND silk sutures. The suture material straightens as it is delivered from the folder. Each suture may be delivered to the surgeon individually from the opening packet or removed from the folder and placed in the suture book. The inner folder for these products is either red with a black C/R symbol or white with red lettering. The safety organizer tray is used for Coated VICRYL sutures, MONOCRYL sutures, PDS II sutures, ETHIBOND EXCEL sutures, PERMA-HAND silk sutures, NUROLON nylon sutures, MERSILENE sutures, and surgical gut sutures. The safety organizer tray allows for single strand arming and dispensing. The needles are situated in individually numbered needle parks and may be armed and dispensed with little or no hand-to-needle contact.

Ligating material used as either single strand (free or freehand) ties, or as continuous ties unwound from a reel or other device. The length of single strand ties is determined by the depth of the wound. In subcutaneous tissue, quarter lengths (approximately 14 inches) are usually long enough for ligating. Deeper wounds may require single strand ties from 18 to 30 inches long. Single strand ligating material is available in pre-cut lengths or 18, 24, and 30 inch strands.

Many surgeons prefer continuous ties. Some prefer LIGAPAK ligature, which is supplied on disc like plastic radiopaque dispensing reels that are color coded by material. The size of the ligature material is indicated by the number of holes visible on the side of the reel (e.g., 3 holes = 3-0 suture). The reel is held in the palm of the hand as blood vessels are ligated. Other surgeons may prefer the ligating material rewound onto a rubber reel, gauze sponge, metal bobbin, or other device.

The number of packets of ligating material required to tie off subcutaneous vessels (bleeders) will vary with patient size and age, the amount of bleeding, the type of operation, the length of the incision, and the surgical technique. An abdominal incision 8 to 12 inches long might require one to three packets to ligate the subcutaneous blood vessels.

All suture material is packaged dry with the exception of surgical gut and pliability ETHILON sutures. Natural absorbable suture materials are packaged with a small amount of sterile fluid, usually alcohol with water, to maintain pliability. They should therefore be opened over a basin to prevent any solution from spilling onto the sterile field.

All needles should be counted after packets of swaged sutures are opened, according to established hospital procedure. The packets should be retained to facilitate verification of the final needle count after the surgical procedure.

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**FIG.7**

### **Suture Handling Technique**

1. With a rolling-outward motion, peel the flaps apart to approximately one-third the way down the sealed edges. Keeping pressure between the knuckles for control, after the

2. Clamp the needleholder approximately one-third to one-half of the distance from the swage area to the needle point. Do not clamp the swaged area. Gently pull the

3. Additional suture straightening should be minimal. If the strand must be straightened, hold the armed needleholder and gently pull the strand making

minimal. If the strand must be straightened, hold the armed needleholder and gently pull the strand making certain not to disarm the needle from the suture.

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## Referencias

1. Artandi C: Industrial sterilization, Point of View 16(2):14-15, April 1979
2. Kohn, ML: Sterile transfer - to flip or not to flip, Point of View 16(2):16, April 1979

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### Hemostasis Products

---

Hemostasis is a complex process that prevents or stops blood flow from vessels in the operative site. Hemostasis not only prevents loss of the patient's blood, but provides as bloodless a field as possible for accurate dissection. Various mechanical, electrical, and topical methods are used to achieve hemostasis.

## Mechanical Method

ETHICON bone wax is a sterile mixture of semisynthetic beeswax and isopropyl palmitate (a wax-softening agent). It is used to control local bleeding from bone surfaces by creating a mechanical (tamponade) barrier.

Bone wax is minimally resorbable and should be used sparingly as it may inhibit osteogenesis and act as a physical barrier to the reparative process. Excess bone wax should be removed from the operative site.

Bone wax should not be used where rapid osseous regeneration and fusion are desired. Mild tissue reaction adjacent to the site of bone wax use has been reported. Also, studies have suggested that bone wax acts as a foreign body and may impair the ability of cancellous bone to clear bacteria.<sup>1</sup> Local accumulation of foreign body giant cells has been observed in animal models. Histologic examination has revealed the appearance of macrophages and occasionally polymorphonuclear leukocytes and lymphocytes.

Bone wax is supplied sterile in 2.5 gram foil envelopes with sealed overwraps. Packets should be opened just prior to use to minimize the possibility of contamination and excessive drying and should be used immediately. Using sterile technique, the wax should be warmed to the desired consistency by manipulation with the fingers or by immersion of the unopened foil packet in warm sterile solution. Bone wax should not be resterilized or subjected to excessive heat.

## Electrical Method

The ETHICON POWERSTAR bipolar scissors provide precise cutting and simultaneous coagulation in open procedures. Hemostasis is achieved through bipolar electrosurgery, or the controlled use of high frequency electric current for cutting and coagulation.

There are two types of electrosurgery: monopolar and bipolar. Both draw high voltage, high frequency current from an electrosurgical generator. The main difference between the two is the method of grounding each uses.

In monopolar electrosurgery, the active electrode is both the cutting and coagulation instrument. The patient becomes part of the return path by being connected to a return electrode in the form of a grounding pad. The patient's tissue provides the impedance in this circuit. However, although current density is diminished by both distance from the contact site and the return electrode, cutting or coagulation current can pass through the patient beyond the area of surgical concern, and may produce thermal damage in tissue surrounding the area of interest. Monopolar electrosurgery lacks a tight focus of operation in patient tissue.

In bipolar electrosurgery, both electrodes are active and are incorporated into a single instrument. Current is focused from one active electrode to the other active electrode through only the very small segment of patient tissue between them. Except for the tissue between the active and return electrodes, the patient is not part of the bipolar circuit. Therefore, bipolar electrosurgery is highly focused with regard to the targeted tissue.

Since electrical energy delivered to tissue by means of bipolar electrosurgery is more focused than is the energy distribution produced by monopolar electrosurgery, cutting and coagulation with the POWERSTAR bipolar scissors can be performed closer to delicate structures, such as nerves and bowel tissue. The bipolar scissors are particularly effective in cardiovascular, peripheral vascular, plastic, and ENT procedures. However, the POWERSTAR bipolar scissors are not indicated for contraceptive coagulation of the fallopian tube, but may be used to achieve hemostasis following transection of the fallopian tube.

POWERSTAR bipolar scissors are available in a variety of styles and sizes and are compatible with most electrosurgical generators. Supplied nonsterile, the reusable POWERSTAR bipolar scissors are sterilized by steam autoclave.

## Topical Agents

---

SURGICEL absorbable hemostat (oxidized regenerated cellulose) is a fabric made from pure alpha cellulose. It is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small artery hemorrhage when ligation or other conventional methods of control are impractical or ineffective. As a result of the regeneration process, oxidized regenerated cellulose absorbs more rapidly than regular oxidized cellulose and is bactericidal.

Upon contact with blood, SURGICEL absorbable hemostat generally stops bleeding within 2 to 8 minutes by acting as a matrix for the formation of a clot. SURGICEL absorbable hemostat is generally absorbed from the site of implantation within 7 to 14 days with a minimum of tissue reaction. It is safe and reliable for use on or near vital tissue. However, SURGICEL hemostat must always be removed from the site of application when used in, around, or in the proximity of bone defects and foramina.

Uniquely, SURGICEL hemostat is bactericidal *in vitro* against a wide range of gram-positive and gram-negative organisms including aerobes and anaerobes. As such, it should not be impregnated with anti-infective agents.

Two variations on regular SURGICEL hemostat, which is a loosely knit fabric, are provided to give the surgeon superior handling qualities for specific applications. SURGICEL NU-KNIT absorbable hemostat is a more densely woven fabric that can be sutured and is especially appropriate for procedures involving solid organs, like the liver. SURGICEL fibrillar absorbable hemostat differs from standard SURGICEL and SURGICEL NU-KNIT in that it has been through a carding/fibrillation process which creates a fibrous web of ORC fibers. This process results in the product having consistency and adherence properties similar to those of cotton.

SURGICEL absorbable hemostat is available in four sizes: 2 x 14 inch (5.1 x 35.6cm), 4 x 8 inch (10.2 x 20.3cm), 2 x 3 inch (5.1 x 7.6cm), and 1/2 x 2 inch (1.3 x 5.1cm). SURGICEL NU-KNIT absorbable hemostat is available in four sizes: 1 x 1 inch (2.5 x 2.5cm), 1 x 3.5 inch (2.5 x 8.9 cm) [endoscopic size], 3 x 4 inch (7.6 x 10.2cm), and 6 x 9 inch (15.2 x 22.9cm). SURGICEL fibrillar is available in three sizes: 1 x 2 inch (2.5 x 5.1cm), 2 x 4 inch (5.1 x 10.2cm), and 4 x 4 inch (10.2 x 10.2cm). All products are supplied sterile and should not be resterilized.

INSTANT collagen absorbable hemostat, a purified and lyophilized (freeze-dried) bovine dermal collagen, is indicated in surgical procedures (other than in neurosurgery and ophthalmology) for use as an adjunct to hemostasis when the control of bleeding by ligature or other conventional methods are ineffective or impractical.

The material, prepared as a spongelike pad, combines the efficacy of collagen to control bleeding with the excellent handling properties of a sponge. INSTANT hemostat, which can be cut to size, is applied directly to the bleeding surface with pressure. Hemostasis is typically achieved in 2 to 5 minutes. INSTANT hemostat maintains its integrity in the presence of blood and is not dispersed when wet. It can be used dry or moistened with sterile saline, depending on surgeon preference.

INSTANT hemostat is supplied sterile in 1 x 2-inch (2.5 x 5.1cm) and 3 x 4-inch (7.6 x 10.2cm) peelable plastic envelopes. It is inactivated by autoclaving and, therefore, should not be resterilized.

INSTANT MCH Microfibrillar Collagen Hemostat, is a dry white absorbent hemostatic agent in microfibrillar form obtained from bovine deep flexor tendon (achilles tendon).

INSTANT MCH (other than neurological, urological, and ophthalmological) is recommended for use in surgical procedures as an adjunct to hemostasis when control of bleeding by ligature or conventional procedures is ineffective or impractical.

The microfibrillar form of the product allows the surgeon to grasp with forceps any amount of INSTANT MCH needed to achieve hemostasis at a particular bleeding site. The microfibrillar form may be more convenient than the sponge form for hard to reach or irregularly shaped bleeding sites. Although it is easy to pull the desired amount of INSTANT MCH from the entire supply, the group of selected fibers continue to cohere to one another and application to the bleeding site is easily controlled.

Topical application of INSTANT MCH effectively controls bleeding, usually within 2 to 4 minutes, when applied directly to the bleeding site. INSTANT MCH is designed to be completely absorbable if left *in situ* after hemostasis. In contact with blood, the fibers expand to become a coherent gelatinous mass that conforms to the shape of the bleeding area. If desired, removal of this mass is easily accomplished.

Sterile INSTAT MCH is supplied in 1 gram and 0.5 gram sizes. It is supplied sterile and should not be resterilized.

THROMBOGEN topical thrombin, U.S.P. is an absorbable hemostatic drug that is used in a wide variety of procedures as an aid in hemostasis whenever oozing blood from capillaries and small venules is accessible.

THROMBOGEN topical thrombin is a protein substance derived from bovine blood and is supplied as a sterile powder that has been lyophilized (freeze dried). The speed with which it will clot the blood is dependent on the concentration of thrombin and diluent. It can be used via spray, splash, powder, absorbable gelatin sponge, or collagen sponge.

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## Tissue Adhesives

Wound closure should be simple, rapid, inexpensive, and painless, and should achieve optimal cosmetic results.<sup>3</sup> The ideal surgical adhesive would be safe for topical application, easy to apply, polymerize rapidly, support the approximated skin edges and maintain the skin edge eversion necessary for maximum wound healing and acceptable cosmesis, and eliminate the need for suture removal.<sup>4</sup>

Tissue adhesives offer many of the advantages of the ideal wound closure devices. They have been shown to be inexpensive and painless, and have antimicrobial activity against gram-positive organisms. They have a low rate of dehiscence and a low infection rate, and provide excellent cosmetic results.<sup>3</sup>

The use of DERMABOND® Topical Skin Adhesive (2-Octyl Cyanoacrylate) significantly decreases the time of treatment for wound closure and eliminates the need for postoperative suture removal. And, although 5- to 7-day follow-up visits may still be necessary, patients will not experience the anxiety and discomfort often associated with suture removal. Additional benefits of DERMABOND adhesive include ease of use, formation of its own protective barrier, eliminating the need for additional bandages, and excellent cosmetic outcome.<sup>4</sup>

Although tissue adhesives have many advantages, successful incorporation into surgical use depends on an understanding of the indications, contraindications, and proper method of application. (*See instructions on next page and Package Insert on pages 118-121.*)

**FIG.2**

### A Guide For Using Dermabond Adhesive


1. Clean and debride the wound and establish hemostasis. Be sure wound edges and surrounding skin are dry, to assure direct tissue contact and prevent premature polymerization of DERMABOND adhesive.

Use subcutaneous sutures as needed. DERMABOND adhesive should not be used as a replacement for proper subcutaneous closure. MONOCRYL® (poliglecaprone 25) suture would be an excellent choice for subcutaneous closure.


2. Remove applicator from packaging and hold with tip pointed upward. Apply pressure at the midpoint of the ampule, crushing the inner glass ampule. Invert applicator and gently squeeze to express the liquid through the applicator tip. After crushing the inner glass ampule, use DERMABOND adhesive immediately.


3. To prevent inadvertent run-off of DERMABOND adhesive, position the wound in a horizontal plane. Manually approximate the wound edges with forceps or gloved fingers. Use gentle brushing strokes to apply a thin film of liquid to the approximated wound edges, and maintain proper eversion of skin edges as you apply DERMABOND adhesive. The adhesive should extend at least 1/2 centimeter on each side of the apposed wound edges. Apply DERMABOND adhesive from above the wound.


Avoid seepage into the wound as it may delay healing.

4. When using DERMABOND adhesive on the face, it is important to prevent the product from trickling into the eyes; failure to do so may seal the eyes shut. Prior to application, apply petroleum jelly around the eye to block the adhesive from entering the eye. Protect and hold the eye closed with a dry gauze pad and position the patient with a slight horizontal tilt so that any run-off travels away from the eye.

5. Gradually build up three or four thin layers of adhesive. Ensure the adhesive is evenly distributed over the wound. Maintain approximation of the wound edges until the adhesive sets and forms a flexible film. This should occur about 1 minute after applying the last layer.

6. Do not apply ointments or medications on top of DERMABOND adhesive.

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## ***Referencias***

- Johnson P, Fromm D: Effects of bone wax on bacterial clearance, *Surgery* 89(2):206-209, 1981
- Barrer S et al: Ideal laparotomy closure: Comparison of retention sutures with new retention bridging devices, *Am Surg* 42:582, 1976
- Singer AJ, et al: Prospective, Randomized, Controlled Trial of Tissue Adesive (2-Octylcyanoacrylate) vs Standard Wound Closure Technique for Laceration Repair, *Acad Emer Med* 5(2):94-99, February 1998
- Toriumi DM, et al: Use of Octyl-2-Cyanoacrylate for Skin Closure in Facial Plastic Surgery, *Plas & Recon Surg* 102(6):2209-2219, November 1998
- Hunt TK, Dunphy JE (ed): *Fundamentals of Wound Management*, New York: Appleton-Century-Crofts, 1979, p. 438
- Myers MB: Sutures and wound healing, *Am J Nurs* 71(9):1726, September 1971
- Stone HH et al: Management of acute full-thickness losses of the abdominal wall, *Ann Surg* 193(5):616, May 1981

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## ***Referencias***

International Standard ISO-10993. Biological Evaluation of Medical Devices. Part 1: Evaluation and Testing.

General Program Memorandum #G95-1, Department of Health and Human Services, Office of Device Evaluation, May 1995.

Guideline for Basic Biological Tests of Medical Materials and Devices. Japan Ministry of Health and Welfare, 1995.

Levy AJ: Testing new suture materials, Point of View, 19(2):15, April 1982

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## Referencias

- Abidin MR et al: Metallurgically bonded needle holder jaws, *Am Surg* 56(10):643-647, October 1990
- American College of Surgeons, Altemeier WA, et al (eds): *Manual on Control of Infection in Surgical Patients*, Philadelphia: Lippincott, 1976
- Amshel AL: The use of Vicryl (polyglactin 910) sutures in colonic and rectal surgery, *Dis Col & Rect* 20(7):635-638, October 1977
- Artandi C: Industrial sterilization, *Point of View* 16(2):14-15, April 1979
- Barrer S et al: Ideal laparotomy closure: Comparison of retention sutures with new retention bridging devices, *Am Surg* 42:582, 1976
- Blaydes JE, Berry J: Comparative evaluation of coated and uncoated polyglactin 910 in cataract and muscle surgery, *Ophth Surg* 11(11):790-793, November 1980
- Chu CC, Williams DF: Effects of physical configuration and chemical structure of suture material on bacterial adhesion, *Am J Surg* 147:197-204, February 1984
- Chusak RB, Dibbell DG: Clinical experience with polydioxanone monofilament absorbable sutures in plastic surgery, *Plas & Recon Surg* 72(2):217-220, August 1983
- Craig PH, et al: A biologic comparison of polyglactin 910 and polyglycolic acid synthetic absorbable sutures, *Surg Gyn & Obs* 14(1):1-10, July 1975
- Cruse P, Foorde R: The epidemiology of wound infection, *Surg Clin NA* 60(1): 27-40, 1980
- Davis JH: *Clinical Surgery*, CV Mosby Co., 1987
- Dineen P: The effect of suture material in the development of vascular infection. *Vasc Surg* 11(1):29-32, January/February 1977
- Durdey P, Bucknall TE: Assessment of sutures for use in colonic surgery: An experimental study, *J Royal Society of Med* 77:472-477, June 1984
- Edlich RF et al: Scientific basis for selecting surgical needles and needle holders for wound closure, *Clin Plastic Surg* 17(3):583-601, July 1990
- Gotttrup F: *Acta Chir Scand* 149:57, 1983
- Herrman JB: Tensile strength and knot security of surgical suture materials, *Am Surg* 37:211, April 1971
- Houdart R et al: Polydioxanone in digestive surgery, *Am J Surg* 152:268-271, September 1986
- Hunt TK, Dunphy JE (ed): *Fundamentals of Wound Management*, New York: Appleton-Century-Crofts, 1979, p. 438
- Hunter P: 10-0 Vicryl suture offers alternative to no stitch, *Ocular Surg News* 9(16): August 15, 1991
- Johnson P, Fromm D: Effects of bone wax on bacterial clearance, *Surgery* 89(2):206-209, 1981
- Kenady DE: Management of abdominal wounds. *Surg Clin NA* 64(4):803-806, 1984
- Kohn ML: Sterile transfer--to flip or not to flip, *Point of View* 16(2):16, April 1979
- Kronenthal RL: Intraocular degradation of nonabsorbable sutures, *Amer Intra-Ocular Implant Soc J* 3:222-228, July/October 1977
- Levy AJ: Testing new suture materials, *Point of View*, 19(2):15, April 1982
- Lo Piccolo M (ed): Mersilene suture results found comparable to nylon, *Ophthalmology Times*:18, August 15, 1988
- McCollum CN, Kester RC: A sliding suture technique for inaccessible arterial anastomoses, *Surg Gynecol Obstet* 153(6):907, 1981
- Murray DH, Blaisdell FW: Use of synthetic absorbable sutures for abdominal and chest wound closure, *Arch Surg* 113:477-480, April 1978
- Myers MB: Sutures and wound healing, *Am J Nurs* 71(9):1726, September 1971
- Nealon TF: *Fundamental Skills in Surgery*, 3d ed, Philadelphia: Saunders, 1979, p. 47
- Pfeffer WH: How to choose sutures for microsurgery, *Contemporary OB/GYN*:1-5, November 1983
- Quinn J, et al: A randomized trial comparing octylcyanoacrylate tissue adhesive and sutures in the management of lacerations, *J Am Med Assn* 277(19):1527-1530, May 21, 1997
- Singer AJ, et al: Prospective, Randomized, Controlled Trial of Tissue Adhesive (2-Octylcyanoacrylate) vs Standard Wound Closure Technique for Laceration Repair, *Acad Emer Med* 5(2):94-99, February 1998
- Stone HH et al: Management of acute full-thickness losses of the abdominal wall, *Ann Surg* 193(5):616, May 1981
- Taylor EJ et al (ed): *Dorland's Medical Dictionary*, Philadelphia: Saunders, 1988

Tera H, Aberg C: Tissue holding power to a single suture in different parts of the alimentary tract, *Acta Chir Scand* 142:343-348, 1976

Thacker JG et al: Surgical needle sharpness, *Am J Surg* 157:334-339, March 1989

Toriumi DM, et al: Use of Octyl-2-Cyanoacrylate for Skin Closure in Facial Plastic Surgery, *Plas & Recon Surg* 102(6):2209-2219, November 1998

Towler MA et al: Laser-drilled cardiovascular surgical needles, *J Appl Biomaterials* 2:183-186, 1991

Trier WC: Considerations in the choice of surgical needles, *Surg Gynecol Obstet* 149:84, 1979

Troutman RC, McGregor W: New compound curved needle TG6-C Plus, *Ophthalmology, Instrument and Book Suppl* 92(8):84-86, 1985

Vallors B, et al: Absorbable or Nonabsorbable Suture Materials for Closure of the Dura Mater?, *Neurosurgery* 9(4):407-413, 1981

VanWinkle Jr W: Epithelization and wound contraction, In *Symposium on Basic Science in Plastic Surgery* (Krizek TJ, Hoopes JE, eds). Mosby, St. Louis 1976, p. 87

VanWinkle Jr W: The healing of skin and visceral wounds and the effect of suture materials, In *Symposium on Basic Science in Plastic Surgery* (Krizek TJ, Hoopes JE, eds). Mosby, St. Louis 1976, p. 135

VanWinkle Jr W, Hastings JC: Considerations in the choice of suture material for various tissues, *Surg Gyn & Obs* 135:113-126, July 1972

# Product Information

## *Surgical Gut Suture*

ABSORBABLE SURGICAL SUTURES, U.S.P.

### **DESCRIPTION**

Surgical gut suture is an absorbable, sterile surgical suture composed of purified connective tissue (mostly collagen) derived from either the serosal layer of beef (bovine) or the submucosal fibrous layer of sheep (ovine) intestines. Surgical gut sutures are available in plain or chromic. Chromic gut is processed to provide greater resistance to absorption. Surgical gut is packaged in tubing fluid. Blue dyed chromic gut suture is also available. Surgical gut suture meets all requirements established by the United States Pharmacopeia (U.S.P.) for absorbable surgical sutures.

### **INDICATIONS**

Surgical gut suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological tissues.

### **ACTIONS**

When surgical gut suture is placed in tissue, a moderate tissue inflammation occurs which is characteristic of foreign body response to a substance. This is followed by a loss of tensile strength and a loss of suture mass, as the proteolytic enzymatic digestive process dissolves the surgical gut. This process continues until the suture is completely absorbed. Many variable factors may affect the rate of absorption. Some of the major factors which can affect tensile strength loss and absorption rates are:

1. Type of suture - plain gut generally absorbs more rapidly than chromic gut.
2. Infection - surgical gut is absorbed more rapidly in infected tissue than in non-infected tissue.
3. Tissue sites - surgical gut will absorb more rapidly in tissue where increased levels of proteolytic enzymes are present, as in the secretions exhibited in the stomach, cervix and vagina.

### **CONTRAINDICATIONS**

This suture, being absorbable, should not be used where extended approximation of tissue is required. The use of this suture is contraindicated in patients with known sensitivities or allergies to collagen or chromium, as gut is a collagen based material, and chromic gut is treated with chromic salt solutions.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving gut suture before using surgical gut suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance when selecting a suture. The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing. As this is an absorbable material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of sites which may undergo expansion, stretching or distention or which may require additional support. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. As an absorbable suture, surgical gut may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. Do not resterilize. Discard open packages and unused sutures. Certain patients may be hypersensitive to collagen or chromium and might exhibit an immunological reaction resulting in inflammation, tissue granulation or fibrosis, wound suppuration and bleeding, as well as sinus formation.

## ***Fast Absorbing Surgical Gut (Plain)***

ABSORBABLE SURGICAL SUTURES, U.S.P.

### **DESCRIPTION**

Fast absorbing surgical gut suture is a strand of collagenous material prepared from the submucosal layers of the small intestine of healthy sheep, or from the serosal layers of the small intestine of healthy cattle. Fast absorbing surgical gut sutures are sterile and elicit only a slight to minimal tissue reaction during absorption. Fast absorbing surgical gut sutures differ from U.S.P. minimum strength requirements by less than 30%.

### **INDICATIONS**

Fast absorbing surgical gut sutures are intended for dermal (skin) suturing only. They should be utilized only for external knot tying procedures.

### **ACTIONS**

The results of implantation studies of fast absorbing surgical gut sutures in the skin of animals indicate that nearly all of its original strength is lost within approximately seven (7) days of implantation. When surgical gut suture is placed in tissue, a moderate tissue inflammation occurs which is characteristic of foreign body response to a substance. This is followed by a loss of tensile strength and a loss of suture mass, as the proteolytic enzymatic digestive process dissolves the surgical gut. This process continues until the suture is completely absorbed. Many variable factors may affect the rate of absorption. Some of the major factors which can affect tensile strength loss and absorption rates are:

1. Type of suture - plain gut generally absorbs more rapidly than chromic gut.
2. Infection - surgical gut is absorbed more rapidly in infected tissue than in non-infected tissue.
3. Tissue sites - surgical gut will absorb more rapidly in tissue where increased levels of proteolytic enzymes are present, as in the secretions exhibited in the stomach, cervix and vagina.

Data obtained from implantation studies in rats show that the absorption of these sutures is essentially complete by the twenty-first (21st) to forty-second (42nd) post implantation day.

### **CONTRAINDICATIONS**

This suture, being absorbable, should not be used where extended approximation of tissue is required. These sutures have been designed to absorb at a rapid rate and must be used on dermal tissue only. These sutures should never be used on internal tissue. The use of this suture is contraindicated in patients with known sensitivities or allergies to collagen, as gut is a collagen based material.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving gut suture before using fast absorbing surgical gut suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing. As this is an absorbable material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of sites which may undergo expansion, stretching or distention or which may require additional support.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. As an absorbable suture, fast absorbing surgical gut may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. Do not resterilize. Discard open packages and unused sutures. Store at room temperature. Certain patients may be hypersensitive to collagen and might exhibit an immunological reaction resulting in inflammation, tissue granulation or fibrosis, wound suppuration and bleeding, as well as sinus formation.

### **PRECAUTIONS**

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Surgical gut sutures require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of

the surgeon. The surgeon should avoid unnecessary tension when running down knots, to reduce the occurrence of surface fraying and weakening of the strand. Avoid prolonged exposure to elevated temperatures. To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

#### **ADVERSE REACTIONS**

Adverse effects associated with the use of this device include wound dehiscence, variable rates of absorption over time (depending on such factors as the type of suture used, the presence of infection and the tissue site), failure to provide adequate wound support in closure of sites where expansion, stretching or distension occur, etc., unless additional support is supplied through the use of nonabsorbable suture material, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from cancer, anemia, obesity, diabetes, infection or other conditions which may delay wound healing, allergic response in patients with known sensitivities to collagen which may result in an immunological reaction resulting in inflammation, tissue granulation or fibrosis, wound suppuration and bleeding, as well as sinus formation, infection, moderate tissue inflammatory response characteristic of foreign body response, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

#### **HOW SUPPLIED**

Fast absorbing surgical gut sutures are available in sizes 5-0 (metric size 1.5) and 6-0 (metric size 1.0) with needles attached in one, two and three dozen boxes.

## ***Coated VICRYL RAPIDE (Polyglactin 910) Braided Coated Synthetic Absorbable Suture, Undyed***

NON-U.S.P.

### **DESCRIPTION**

Coated VICRYL RAPIDE (polyglactin 910) suture is a synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide. The empirical formula of the copolymer is  $(C_2H_2O_2)_m(C_3H_4O_2)_n$ . The characteristic of rapid loss of strength is achieved by use of a polymer material with a lower molecular weight than Coated VICRYL (polyglactin 910) suture. Coated VICRYL RAPIDE sutures are obtained by coating the braided suture material with a mixture of equal parts of copolymer of glycolide and lactide (polyglactin 370) and calcium stearate. Polyglactin 910 copolymer and polyglactin 370 with calcium stearate have been found to be nonantigenic, nonpyrogenic and elicit only a mild tissue reaction during absorption. Coated VICRYL RAPIDE sutures are only available undyed. Although this suture is a synthetic absorbable suture, its performance characteristics are intended to model the performance of collagen (surgical gut) suture. The knot tensile strength of Coated VICRYL RAPIDE suture meets U.S.P. knot tensile strength requirements for collagen sutures, however, Coated VICRYL RAPIDE suture strength is up to 26% less than knot tensile strength requirements for synthetic absorbable sutures.

---

### **MAXIMUM SUTURE OVERSIZE IN DIAMETER (MM) FROM U.S.P**

---

U.S.P SUTURE SIZE DISIGNATION	MAXIMUM OVERSIZE (MM)
5-0	.016
4-0	.017
3-0	.018
2-0	.010
0	.022

---

### **INDICATIONS**

Coated VICRYL RAPIDE synthetic absorbable suture is indicated only for use in superficial soft tissue approximation of the skin and mucosa, where only short term wound support (7-10 days) is required. Coated VICRYL RAPIDE suture is not intended for use in ligation, ophthalmic, cardiovascular or neurological procedures.

### **ACTIONS**

Coated VICRYL RAPIDE suture, when used in closure of skin and mucous membranes, typically begins to fall off 7-10 days post-operatively and can be wiped off subsequently with sterile gauze. Natural mechanical abrasion of the sutures while in situ may also accelerate this disappearance rate. Rapid loss of tensile strength may preclude the need for stitch removal. Coated VICRYL RAPIDE elicits a minimal to moderate acute inflammatory reaction in tissue. Progressive loss of tensile strength and eventual absorption of Coated VICRYL RAPIDE occurs by means of hydrolysis, where the copolymer degrades to glycolic and lactic acids which are subsequently absorbed and metabolized in the body. Absorption begins as a loss of tensile strength followed by a loss of mass. Subcutaneous tissue implantation studies of Coated VICRYL RAPIDE sutures in rats show that 5 days post-implantation approximately 50% of the original tensile strength remains. All of the original tensile strength is lost by approximately 10 to 14 days post-implantation. Intramuscular implantation studies in rats show that the absorption of these sutures occurs thereafter and is essentially complete by 42 days.

### **CONTRAINDICATIONS**

Due to the rapid loss of tensile strength, this suture should not be used where extended approximation of tissues under stress is required or where wound support beyond 7 days is required.

## **WARNINGS**

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing Coated VICRYL RAPIDE suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance when selecting a suture. The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing. Do not resterilize. Discard opened packages and unused sutures. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. As an absorbable suture, Coated VICRYL RAPIDE suture may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. As this is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention, or which may require additional support.

## **PRECAUTIONS**

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Coated VICRYL RAPIDE suture, which is treated with coating to enhance handling characteristics, requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. Avoid prolonged exposure to elevated temperatures. To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

## **ADVERSE REACTIONS**

Adverse effects associated with the use of this device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

## **HOW SUPPLIED**

Coated VICRYL RAPIDE sutures are available sterile, undyed and attached to stainless steel needles of varying types and sizes. Coated VICRYL RAPIDE sutures are available in various lengths in sizes 5-0 to 1 (1.0 to 4.0 metric) in one and three dozen boxes.

## **Coated VICRYL (Polyglactin 910) Suture**

U.S.P., EXCEPT FOR DIAMETER DESCRIPTION Coated VICRYL\* (polyglactin 910) suture is a synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide. Coated VICRYL suture is prepared by coating VICRYL suture material with a mixture composed of equal parts of copolymer of glycolide and lactide (polyglactin 370) with calcium stearate have been found to be nonantigenic, nonpyrogenic and elicits only a mild tissue reaction during absorption. The sutures are available dyed and undyed (natural). Coated VICRYL sutures are U.S.P. except for diameters in the following sizes:

---

### **MAXIMUM SUTURE OVERSIZE IN DIAMETER (MM) FROM U.S.P**

---

U.S.P SUTURE SIZE DISIGNATION	MAXIMUM OVERSIZE (MM)
6-0	.008
5-0	.016
4-0	.017
3-0	.018
2-0	.010
0	.022

---

### **INDICATIONS**

Coated VICRYL suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular or neurological tissues.

### **ACTIONS**

Coated VICRYL suture elicits a minimal acute inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of Coated VICRYL suture occurs by means of hydrolysis, where the copolymer degrades to glycolic and lactic acids which are subsequently absorbed and metabolized in the body. Absorption begins as a loss of tensile strength followed by a loss of mass. Implantation studies in rats indicate that Coated VICRYL suture retains approximately 75% of the original tensile strength at two weeks post implantation. At three weeks, approximately 50% of the original strength is retained for sizes 6-0 and larger and approximately 40% of its original strength is retained for sizes 7-0 and smaller. All of the original tensile strength is lost between four and five weeks post implantation. Absorption of Coated VICRYL suture is essentially complete between 56 and 70 days.

---

Days Implantation	Approximate % Original Strength Remaining
14 Days	75%
21 Days (6-0 and larger)	50%
21 Days (7-0 and Smaller)	40%

---

### **CONTRAINDICATIONS**

This suture, being absorbable, should not be used where extended approximation of tissue under stress is required.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing Coated VICRYL suture for wound closure, as risk of wound dehiscence may vary with the site

of application and the suture material used. Physicians should consider the in vivo performance (under ACTIONS section) when selecting a suture.

The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing. As this is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention, or which may require additional support.

Do not resterilize. Discard opened packages and unused sutures.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Coated VICRYL suture may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

#### **PRECAUTIONS**

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated.

Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon.

Consideration should be taken in the use of absorbable sutures in tissue with poor blood supply as suture extrusion and delayed absorption may occur.

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Coated VICRYL sutures, which are treated to enhance handling characteristics, require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon.

Avoid prolonged exposure to elevated temperatures.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third ( $\frac{1}{3}$ ) to one-half ( $\frac{1}{2}$ ) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

#### **ADVERSE REACTIONS**

Adverse effects associated with the use of this device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

#### **HOW SUPPLIED**

Coated VICRYL sutures are available sterile, as braided dyed (violet) and undyed (natural) strands in sizes 8-0 through 3 (metric sizes 0.4-

## **MONOCRYL Violet Monofilament (Poliglecaprone 25) Suture**

SYNTHETIC ABSORBABLE SUTURE, U.S.P., EXCEPT FOR DIAMETER

### **DESCRIPTION**

MONOCRYL (poliglecaprone 25) suture is a monofilament synthetic absorbable surgical suture prepared from a copolymer of glycolide and epsilon-caprolactone. Poliglecaprone 25 copolymer has been found to be nonantigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption.

MONOCRYL sutures are U.S.P. except for diameters in the following sizes:

---

### **MAXIMUM SUTURE OVERSIZE IN DIAMETER (MM) FROM U.S.P.**

---

U.S.P. SUTURE SIZE DESIGNATION	MAXIMUM OVERSIZE (MM)
6-0	0.049
5-0	0.033
4-0	0.045
3-0	0.067
2-0	0.055
0	0.088
1	0.066
2	0.099

---

### **INDICATIONS**

MONOCRYL sutures are indicated for use in general soft tissue approximation and/or ligation, but not for use in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery.

### **ACTIONS**

MONOCRYL suture is a monofilament which elicits a minimal acute inflammatory reaction in tissues and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of MONOCRYL sutures occurs by means of hydrolysis. Absorption begins as a loss of tensile strength followed by a loss of mass. Implantation studies in rats indicate that MONOCRYL suture retains approximately 60 to 70% of its original tensile strength 7 days post implantation, and approximately 30 to 40% of its original tensile strength at 14 days post implantation. Essentially all of the original tensile strength is lost by 28 days post implantation. Absorption of MONOCRYL absorbable synthetic suture is essentially complete between 91 and 119 days.

---

DAYS IMPLANTATION	APPROXIMATE % ORIGINAL STRENGTH REMAINING
7 DAYS	60 TO 70%
14 DAYS	30 TO 40%

---

### **CONTRAINDICATIONS**

This suture, being absorbable, should not be used where extended approximation of tissue under stress is required.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing MONOCRYL suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under ACTIONS section) when selecting a suture for use in patients. The use of this suture may be inappropriate

in elderly, malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing.

Do not resterilize. Discard open packages and unused sutures. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. As an absorbable suture, MONOCRYL suture may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention, or which may require additional support.

#### **PRECAUTIONS**

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with absorption.

Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon.

Consideration should be taken in the use of absorbable sutures in tissue with poor blood supply as suture extrusion and delayed absorption may occur.

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. MONOCRYL suture knots must be properly placed to be secure. Adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

Avoid prolonged exposure to elevated temperatures.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

#### **ADVERSE REACTIONS**

Adverse effects associated with the use of synthetic absorbable sutures include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

#### **HOW SUPPLIED**

MONOCRYL sutures are available as sterile, monofilament, dyed (violet) strands in sizes 6-0 through 2 (metric sizes 0.7-5), in a variety of lengths, with or without needles. MONOCRYL sutures are also available in sizes 3-0 through 1 (metric sizes 2-4) attached to CONTROL RELEASE removable needles. MONOCRYL sutures are available in one and three dozen boxes.

## ***PDS II (Polydioxanone) Suture Dyed and Clear Monofilament***

SYNTHETIC ABSORBABLE SUTURES, U.S.P., EXCEPT FOR DIAMETER

### **DESCRIPTION**

PDS II (polydioxanone) monofilament synthetic absorbable suture is prepared from the polyester, poly (p-dioxanone). The empirical molecular formula of the polymer is (C<sub>4</sub>H<sub>6</sub>O<sub>3</sub>)<sub>x</sub>.

Polydioxanone polymer has been found to be nonantigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption.

PDS II sutures are U.S.P., except for diameter.

---

### **MAXIMUM SUTURE OVERSIZE IN DIAMETER (MM) FROM U.S.P.**

---

U.S.P. SUTURE SIZE DESIGNATION	MAXIMUM OVERSIZE (MM)
9-0	.005
8-0	.008
7-0	.020
6-0	.015
5-0	.029
4-0	.029
3-0	.056
2-0	.029
0	.071
1	.047
2	.023

---

### **ACTIONS**

Two important characteristics describe the in vivo performance of absorbable sutures: first, tensile strength retention, and second, the absorption rate (loss of mass). PDS II synthetic absorbable suture has been formulated to minimize the variability of these characteristics and to provide wound support through an extended healing period.

The results of implantation studies of PDS II monofilament suture in animals indicate that approximately 70% of its original strength remains two weeks after implantation. At four weeks post-implantation, approximately 50% of its original strength is retained, and at six weeks, approximately 25% of the original strength is retained.

Data obtained from implantation studies in rats show that the absorption of these sutures is minimal until about the 90th post-implantation day. Absorption is essentially complete within six months.

### **INDICATIONS**

PDS II monofilament synthetic absorbable sutures are indicated for use in all types of soft tissue approximation, including use in pediatric cardiovascular tissue where growth is expected to occur and ophthalmic surgery. PDS II suture is not indicated in adult cardiovascular tissue, microsurgery and neural tissue. These sutures are particularly useful where the combination of an absorbable suture and extended wound support (up to six weeks) is desirable.

### **CONTRAINDICATIONS**

These sutures, being absorbable, are not to be used where prolonged (beyond six weeks) approximation of tissues under stress is required and are not to be used in conjunction with prosthetic devices, i.e., heart valves or synthetic grafts.

### **WARNINGS**

The safety and effectiveness of PDS II (polydioxanone) sutures have not been established in neural tissue, adult cardiovascular tissue or for use in microsurgery.

Under certain circumstances, notably orthopaedic procedures, immobilization by external support may be employed at the discretion of the surgeon.

Do not resterilize.

### **PRECAUTIONS**

The PDS II suture knots must be properly placed to be secure. As with other synthetic sutures, knot security requires the standard surgical technique of flat and square ties with additional throws if indicated by surgical circumstance and the experience of the operator.

As with any suture, care should be taken to avoid damage when handling. Avoid the crushing or crimping application of surgical instruments, such as needle holders and forceps, to the strand except when grasping the free end of the suture during an instrument tie.

Conjunctival and vaginal mucosal sutures remaining in place for extended periods may be associated with localized irritation and should be removed as indicated.

Subcuticular sutures should be placed as deeply as possible in order to minimize the erythema and induration normally associated with absorption.

Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds.

### **ADVERSE REACTIONS**

Due to prolonged suture absorption, some irritation and bleeding has been observed in the conjunctiva and mild irritation has been observed in the vaginal mucosa.

### **DOSAGE AND ADMINISTRATION**

Use as required per surgical procedure.

### **HOW SUPPLIED**

PDS II sutures are available as sterile, monofilament dyed (violet) strands in sizes 9-0 thru 2 (metric sizes 0.3-5), and sterile, monofilament dyed (blue) strands in sizes 9-0 thru 7-0 (metric sizes 0.3-0.5) in a variety of lengths, with a variety of needles.

PDS II monofilament dyed (violet) sutures, sizes 4-0 thru 1 (metric sizes 1.5-4) are also available attached to CONTROL RELEASE removable needles.

PDS II Clear suture strands are available in sizes 7-0 thru 1 (metric size 0.5-4) in a variety of lengths with permanently attached needles.

CAUTION: Federal Law restricts this device to sale, distribution and use by or on the order of a physician or a veterinarian.

## **PANACRYL Suture**

SYNTHETIC ABSORBABLE SUTURE, U.S.P.,  
EXCEPT FOR DIAMETER

### **DESCRIPTION**

PANACRYL suture is a synthetic, braided, undyed, absorbable, sterile surgical suture composed of a copolymer of 95% lactide and 5% glycolide. The suture is coated with a copolymer composed of 90% caprolactone and 10% glycolide. The copolymers in this product have been shown to be nonantigenic, nonpyrogenic and to elicit mild tissue reaction during absorption.

PANACRYL sutures are U.S.P. except for diameters in the following sizes:

---

### **MAXIMUM SUTURE OVERSIZE IN DIAMETER (MM) FROM U.S.P.**

---

<b>U.S.P. SUTURE SIZE DESIGNATION</b>	<b>MAXIMUM OVERSIZE (MM)</b>
3-0	0.067
2-0	0.055
0	0.088
1	0.066
2	0.099

---

### **INDICATIONS**

PANACRYL suture is indicated for use in general soft tissue approximation and/or ligation, and orthopaedic uses including tendon and ligament repairs and reattachment to bone but not for use in ophthalmic, cardiovascular or neurological tissue. PANACRYL suture is particularly useful where extended wound support (up to 6 months) is desirable.

### **ACTIONS**

PANACRYL suture elicits a minimal inflammatory reaction in tissue and allows ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of PANACRYL suture occurs by means of hydrolysis, where the copolymer degrades to glycolic and lactic acids which are subsequently absorbed and metabolized in the body. Absorption begins as a loss of tensile strength followed by a loss of mass.

Implantation studies in rats indicate that PANACRYL suture retains approximately 80% of its original strength at 3 months and 60% of its original strength at 6 months post implantation. PANACRYL suture is essentially absorbed in 1.5 years or longer.

---

<b>DAYS IMPLANTATION</b>	<b>APPROXIMATE % ORIGINAL STRENGTH REMAINING</b>
3 months	80%
6 months	60%

---

### **CONTRAINDICATIONS**

PANACRYL suture, being absorbable, should not be used where extended approximation of tissue beyond 6 months is required.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing PANACRYL suture for wound closure, as a risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under

ACTIONS Section) when selecting the suture. Use of this suture may be inappropriate in patients suffering conditions which could delay wound healing beyond 6 months.

Do not resterilize. Discard opened packages and unused sutures.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. As an absorbable suture, PANACRYL suture may act transiently as a foreign body.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

#### **PRECAUTIONS**

Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon.

Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur.

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. PANACRYL sutures, which are treated to enhance handling characteristics, require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and the experience of the surgeon.

To avoid damaging needle points and swage areas, grasp the needles in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

Avoid prolonged exposure to elevated temperatures.

#### **ADVERSE REACTIONS**

Adverse effects associated with the use of PANACRYL suture include wound dehiscence, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than seven days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine and bile occurs and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

#### **HOW SUPPLIED**

PANACRYL sutures are available as sterile, braided, undyed (white) strands in sizes 3-0 through 2 in a variety of lengths with and without needles. PANACRYL sutures are also available attached to CONTROL RELEASE removable needles.

PANACRYL sutures are available in one, two, or three dozen boxes.

## **PERMA-HAND Silk Suture**

NONABSORBABLE SURGICAL SUTURE, U.S.P.

### **DESCRIPTION**

PERMA-HAND silk suture is a nonabsorbable, sterile, surgical suture composed of an organic protein called fibroin. This protein is derived from the domesticated species Bombyx mori (B. More) of the family Bombycidae. PERMA-HAND sutures are processed to remove the natural waxes and gums. PERMA-HAND suture is dyed black and coated with a special wax mixture. PERMA-HAND suture is also available in its natural color. PERMA-HAND Virgin silk suture is available in which the sericin gum is not removed and serves to hold the filaments together.

PERMA-HAND suture meets requirements established by the United States Pharmacopoeia (U.S.P.) for nonabsorbable surgical suture.

### **INDICATIONS**

PERMA-HAND suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

### **ACTIONS**

PERMA-HAND suture elicits an acute inflammatory reaction in tissue, which is followed by a gradual encapsulation of the suture by fibrous connective tissue. While silk sutures are not absorbed, progressive degradation of the proteinaceous silk fiber in vivo may result in gradual loss of all of the suture's tensile strength over time.

### **CONTRAINDICATIONS**

The use of this suture is contraindicated in patients with known sensitivities or allergies to silk.

Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, silk should not be used where permanent retention of tensile strength is required.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing PERMA-HAND suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Acceptable surgical practice should be followed for the management of infected or contaminated wounds.

Do not resterilize. Discard opened packages and unused sutures.

### **PRECAUTIONS**

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

As with any suture material, adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

### **ADVERSE REACTIONS**

Adverse effects associated with the use of this device include wound dehiscence, gradual loss of all tensile strength over time, allergic response in patients that are known to be sensitive to silk, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, acute inflammatory tissue reaction, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

### **HOW SUPPLIED**

PERMA-HAND sutures are available in U.S.P. sizes 9-0 through 5 (metric sizes 0.3-7.0) in a variety of lengths with and without permanently attached needles and on LIGAPAK dispensing reels.

PERMA-HAND sutures are also available in U.S.P. sizes 4-0 through 1 (metric sizes 1.5-4.0) attached to CONTROL RELEASE removable needles.

PERMA-HAND sutures are available in one, two, and three dozen boxes.

## ***Surgical Stainless Steel Suture***

NONABSORBABLE SURGICAL SUTURES, U.S.P.

### **DESCRIPTION**

Surgical stainless steel suture is a nonabsorbable, sterile surgical suture composed of 316L stainless steel. Surgical stainless steel suture is available as a monofilament and multifilament suture.

Surgical stainless steel suture meets all requirements established by the United States Pharmacopoeia (U.S.P.) for nonabsorbable, surgical sutures. Surgical stainless steel suture is also labeled with the B&S gauge classifications.

### **INDICATIONS**

Surgical stainless steel suture is indicated for use in abdominal wound closure, hernia repair, sternal closure and orthopaedic procedures including cerclage and tendon repair.

### **ACTIONS**

Surgical stainless steel suture elicits a minimal acute inflammatory reaction in tissue and is not absorbed.

### **CONTRAINDICATIONS**

The use of this suture is contraindicated in patients with known sensitivities or allergies to 316L stainless steel, or constituent metals such as chromium and nickel.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving nonabsorbable, stainless steel sutures before employing for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Acceptable surgical practice must be followed for the management of contaminated or infected wounds.

### **PRECAUTIONS**

In handling this or any other suture material, care should be taken to avoid damage from handling, such as kinking or excessive twisting.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

### **ADVERSE REACTIONS**

Adverse effects associated with the use of this device include wound dehiscence, gradual loss of all tensile strength over time, allergic response in patients with known sensitivities to 316L stainless steel, or constituent metals such as chromium and nickel, infection, minimal acute inflammatory tissue reaction, pain, edema and local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

### **HOW SUPPLIED**

Surgical stainless steel sutures are available in sizes 7 through 10-0 (metric sizes 9.0-0.2) in a variety of lengths with and without permanently attached needles in one, two and three dozen boxes.

## ***ETHILON Nylon Suture***

NONABSORBABLE SURGICAL SUTURES, U.S.P.

### **DESCRIPTION**

ETHILON nylon suture is a nonabsorbable, sterile surgical monofilament suture composed of the long-chain aliphatic polymers Nylon 6 and Nylon 6,6. ETHILON sutures are dyed black or green to enhance visibility in tissue. The suture is also available undyed (clear).

ETHILON suture meets all requirements established by the United States Pharmacopoeia (U.S.P.) for nonabsorbable surgical suture.

### **INDICATIONS**

ETHILON suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

### **ACTIONS**

ETHILON suture elicits a minimal acute inflammatory reaction in tissue, which is followed by gradual encapsulation of the suture by fibrous connective tissue. While nylon is not absorbed, progressive hydrolysis of the nylon in vivo may result in gradual loss over time of tensile strength.

### **CONTRAINDICATIONS**

Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, nylon suture should not be used where permanent retention of tensile strength is required.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing ETHILON suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

Do not resterilize. Discard open packages and unused sutures.

### **PRECAUTIONS**

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

As with any suture material, adequate knot security requires the accepted surgical technique of flat and square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

### **ADVERSE REACTIONS**

Adverse effects associated with the use of this device include wound dehiscence, gradual loss of tensile strength over time, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal acute inflammatory tissue reaction, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

### **HOW SUPPLIED**

ETHILON sutures are available as sterile monofilament strands in U.S.P. sizes 11-0 through 2 (metric sizes 0.1-5.0) in a variety of lengths, with and without permanently attached needles. ETHILON sutures are available in one, two and three dozen boxes.

## **NUROLON Nylon Suture**

NONABSORBABLE SURGICAL SUTURE, U.S.P.

### **DESCRIPTION**

NUROLON nylon suture is a nonabsorbable sterile surgical braided suture composed of the long-chain aliphatic polymers Nylon 6 or Nylon 6,6. NUROLON sutures are dyed black to enhance visibility in tissue. The suture is also available undyed (clear).

NUROLON suture meets all requirements established by the United States Pharmacopoeia (U.S.P.) for nonabsorbable surgical suture.

### **INDICATIONS**

NUROLON suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

### **ACTIONS**

NUROLON suture elicits a minimal acute inflammatory reaction in tissue, which is followed by a gradual encapsulation of the suture by fibrous connective tissue. While nylon is not absorbed, progressive hydrolysis of the nylon in vivo may result in gradual loss of tensile strength over time.

### **CONTRAINDICATIONS**

Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, nylon suture should not be used where permanent retention of tensile strength is required.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing NUROLON suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Acceptable surgical practice should be followed for the management of infected or contaminated wounds.

Do not resterilize. Discard opened packages and unused sutures.

### **PRECAUTIONS**

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical technique of flat and square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

### **ADVERSE REACTIONS**

Adverse effects associated with the use of this device include wound dehiscence, gradual loss of tensile strength over time, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal acute inflammatory tissue reaction, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

### **HOW SUPPLIED**

NUROLON sutures are available in U.S.P. sizes 6-0 through 1 (metric sizes 0.7-4.0) in a variety of lengths with and without permanently attached needles.

NUROLON sutures are available in U.S.P. sizes 4-0 through 1 (metric sizes 1.5-4.0) attached to CONTROL RELEASE removable needles.

NUROLON sutures are available in one, two and three dozen boxes.

## **MERSILENE Polyester Fiber Suture**

NONABSORBABLE SURGICAL SUTURE, U.S.P.

EXCEPT FOR SIZE 6-0 DIAMETER

### **DESCRIPTION**

MERSILENE polyester suture is a nonabsorbable, braided, sterile, surgical suture composed of Poly (ethylene terephthalate). It is prepared from fibers of high molecular weight, long-chain, linear polyesters having recurrent aromatic rings as an integral component. MERSILENE sutures are braided for optimal handling properties, and for good visibility in the surgical field, are dyed green.

Size 6-0 MERSILENE sutures are U.S.P., except for diameter.

---

### **MAXIMUM SUTURE OVERSIZE IN DIAMETER (MM) FROM U.S.P.**

---

U.S.P. SUTURE SIZE DESIGNATION	MAXIMUM OVERSIZE (MM)
6-0	0.024

---

### **INDICATIONS**

MERSILENE suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

### **ACTIONS**

MERSILENE suture elicits a minimal acute inflammatory reaction in tissue, followed by a gradual encapsulation of the suture by fibrous connective tissue. Implantation studies in animals show no meaningful decline in polyester suture strength over time. The polyester fiber suture material is pharmacologically inactive.

### **CONTRAINDICATIONS**

None known.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing MERSILENE suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Do not resterilize. Discard opened packages and unused sutures.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Acceptable surgical practice should be followed for the management of infected or contaminated wounds

### **PRECAUTIONS**

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. The use of additional throws is particularly appropriate when knotting monofilament sutures.

As with any suture material, adequate knot security requires the accepted surgical technique of flat and square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

### **ADVERSE REACTIONS**

Adverse effects associated with the use of this device include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal acute inflammatory tissue reaction and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

### **HOW SUPPLIED**

MERSILENE sutures are available as sterile, braided, green and undyed (white) strands in sizes 6-0 through 5 (metric sizes 0.7-7) in a variety of lengths, with and without permanently attached needles. MERSILENE sutures are also available in green monofilament in U.S.P. sizes 10-0 and 11-0 (metric sizes 0.2-0.1).

MERSILENE sutures, green, braided in U.S.P. size 0 (metric size 3.5) are also available attached to CONTROL RELEASE removable needles.

MERSILENE sutures are available in one, two and three dozen boxes.

## **ETHIBOND Excel Polyester Suture**

NONABSORBABLE SURGICAL SUTURE, U.S.P. EXCEPT FOR SIZE 6-0 DIAMETER

### **DESCRIPTION**

ETHIBOND *EXCEL* polyester suture is a nonabsorbable, braided, sterile, surgical suture composed of Poly (ethylene terephthalate). It is prepared from fibers of high molecular weight, long-chain, linear polyesters having recurrent aromatic rings as an integral component. ETHIBOND *EXCEL* suture is uniformly coated with polybutylate or poly {1xy-1, 4 butanediyl oxy (1, 6-dioxo-1, 6 hexanediyl)}. The highly adherent coating is a relatively nonreactive nonabsorbable compound which acts as a lubricant to mechanically improve the physical properties of the uncoated suture by improving handling qualities as contrasted to the braided, uncoated fiber.

ETHIBOND *EXCEL* sutures are braided for optimal handling properties, and for good visibility in the surgical field, are dyed green.

Size 6-0 ETHIBOND *EXCEL* sutures are U.S.P., except for diameter.

---

### **MAXIMUM SUTURE OVERSIZE IN DIAMETER (MM) FROM U.S.P.**

---

<b>U.S.P. SUTURE SIZE DESIGNATION</b>	<b>MAXIMUM OVERSIZE (MM)</b>
6-0	.0024

---

### **INDICATIONS**

ETHIBOND *EXCEL* suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

### **ACTIONS**

ETHIBOND *EXCEL* suture elicits a minimal acute inflammatory reaction in tissue, followed by a gradual encapsulation of the suture by fibrous connective tissue. Implantation studies in animals show no meaningful decline in polyester suture strength over time. Both polyester fiber suture material and the polybutylate coating are pharmacologically inactive.

### **CONTRAINDICATIONS**

None known.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing ETHIBOND *EXCEL* suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Do not resterilize. Discard opened packages and unused sutures.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Acceptable surgical practice should be followed for the management of infected or contaminated wounds.

### **PRECAUTIONS**

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

As with any suture material, adequate knot security requires the accepted surgical technique of flat and square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

### **ADVERSE REACTIONS**

Adverse effects associated with the use of this device include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs,

infection, minimal acute inflammatory tissue reaction and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

**HOW SUPPLIED**

ETHIBOND *EXCEL* sutures are available as sterile, braided, green and undyed (white) strands in sizes 7-0 through 5 (metric sizes 0.5-7) in a variety of lengths, with and without permanently attached needles.

ETHIBOND *EXCEL* sutures, green, braided, in sizes 4-0 through 1 (metric sizes 1.5-4) are also available attached to CONTROL RELEASE\* removable needles.

ETHIBOND *EXCEL* sutures, green and undyed, are also available attached to TFE polymer pledgets measuring 1/8" x 1/8" x 1/16" (3.0mm x 3.0mm x 1.5mm), 1/4" x 1/8" x 1/16" (7.0mm x 3.0mm x 1.5mm).

ETHIBOND *EXCEL* sutures are available in one, two and three dozen boxes.

## ***PROLENE Polypropylene Suture***

NONABSORBABLE SURGICAL SUTURE, U.S.P.

Except for size 7-0 diameter and HEMO-SEAL Needle Suture Attachment

### **DESCRIPTION**

PROLENE polypropylene suture (clear or pigmented) is a nonabsorbable sterile surgical suture composed of an isotactic crystalline stereoisomer of polypropylene, a synthetic linear polyolefin. The suture is pigmented blue to enhance visibility.

Size 7-0 PROLENE sutures are U.S.P., except for diameter.

---

### **MAXIMUM SUTURE OVERSIZE IN DIAMETER (MM) FROM U.S.P.**

---

<b>U.S.P. SUTURE SIZE DESIGNATION</b>	<b>MAXIMUM OVERSIZE (MM)</b>
7-0	.007

---

PROLENE suture, available as HEMO-SEAL needle suture, is a needle suture combination in which the diameter of the needle swage area has been reduced to facilitate attachment of finer wire diameter needles. The diameter of the suture strand and the needle wire have been more closely aligned to reduce the degree of needle hole bleeding. HEMO-SEAL needle suture differs from U.S.P. in needle attachment requirements only.

### **INDICATIONS**

PROLENE suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

### **ACTIONS**

PROLENE suture elicits a minimal acute inflammatory reaction in tissue, which is followed by a gradual encapsulation of the suture by fibrous connective tissue. PROLENE suture is not absorbed, nor is it subject to degradation or weakening by the action of tissue enzymes. As a monofilament, PROLENE suture, U.S.P. resists involvement in infection and has been successfully employed in contaminated and infected wounds to eliminate or minimize later sinus formation and suture extrusion. The lack of adherence to tissues has facilitated the use of PROLENE suture as a pull-out suture.

### **CONTRAINDICATIONS**

None known.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing PROLENE suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Do not resterilize. Discard opened packages and unused sutures.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Acceptable surgical practice must be followed for the management of infected or contaminated wounds.

### **PRECAUTIONS**

In handling this suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Adequate knot security requires the accepted surgical technique of flat, square ties of single suture strands. The use of additional throws is particularly appropriate when knotting polypropylene sutures.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

## ADVERSE REACTIONS

Adverse effects associated with the use of this device include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal acute inflammatory tissue reaction, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

## HOW SUPPLIED

PROLENE sutures, pigmented, are available as sterile strands in U.S.P. sizes 10-0 through 8-0 (metric sizes 0.2-0.4) and 6-0 through 2 (metric sizes 0.7-5.0).

PROLENE sutures, clear, are available as sterile strands in U.S.P. sizes 6-0 through 2 (metric sizes 0.7-5.0). Size 7-0 (metric size 0.5) PROLENE sutures, pigmented and clear are U.S.P. except for diameter. All PROLENE sutures are available in a variety of lengths, with permanently attached needles.

PROLENE sutures, pigmented and clear are also available as sterile strands in U.S.P. sizes 0 through 2 (metric sizes 3.5-5.0) attached to CONTROL RELEASE removable needles.

PROLENE sutures, pigmented and clear are also available as sterile strands in U.S.P. sizes 0 through 5-0, attached to TFE pledgets measuring 1/4" x 1/8" x 1/16" (7.0mm x 3.0mm x 1.5mm).

PROLENE sutures, pigmented and clear are also available in sterile strands as HEMO-SEAL needle sutures in the following sizes:

PROLENE SUTURE U.S.P. SIZE	HEMO-SEAL NEEDLE SUTURE	HEMO-SEAL NEEDLE SUTURE LIMITS ON NEEDLE ATTACHMENT
	AVG. (KGF)	INDIVIDUAL (KGF) MIN.
5-0	HS-7 0.17	0.08
4-0	HS-6 0.23	0.11
3-0	HS-5 0.45	0.23
USP LIMITS ON NEEDLE ATTACHMENT		
U.S.P. SIZE	AVG. (KGF) MIN.	INDIVIDIAL (KGF) MIN.
5-0	0.23	0.11
4-0	0.45	0.23
3-0	0.68	0.34

PROLENE sutures are available in one, two and three dozen boxes.

## ***PRONOVA Poly (Hexafluoropropylene-VDF) Suture***

NONABSORBABLE SURGICAL SUTURE, U.S.P.

EXCEPT FOR SIZE 7-0 DIAMETER

### **DESCRIPTION**

PRONOVA suture (clear or pigmented) is a nonabsorbable, sterile surgical suture made from a polymer blend of poly(vinylidene fluoride) and poly(vinylidene fluoride-cohexafluoropropylene). The suture is pigmented blue to enhance visibility.

Size 7-0 PRONOVA sutures are U.S.P., except for diameter.

---

### **MAXIMUM SUTURE OVERSIZE IN DIAMETER (MM) FROM U.S.P.**

---

U.S.P. SUTURE SIZE DESIGNATION	MAXIMUM OVERSIZE (MM)
7-0	.007

---

### **INDICATIONS**

PRONOVA suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

### **ACTIONS**

PRONOVA suture elicits a minimal to mild inflammatory reaction in tissue, which is followed by a gradual encapsulation of the suture by fibrous connective tissue. PRONOVA suture is not absorbed, nor is it subject to degradation or weakening by the action of tissue enzymes. As a monofilament, PRONOVA suture, U.S.P. resists involvement in infection and has been successfully employed in contaminated and infected wounds to eliminate or minimize later sinus formation and suture extrusion. Furthermore, the lack of adherence to tissues has facilitated the use of PRONOVA suture as a pull-out suture.

### **CONTRAINDICATIONS**

None known.

### **WARNINGS**

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing PRONOVA suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Do not resterilize. Discard opened packages and unused sutures.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Acceptable surgical practice must be followed for the management of infected or contaminated wounds.

### **PRECAUTIONS**

In handling this suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Adequate knot security requires the accepted surgical technique of flat, square ties of single suture strands.

The use of additional throws is particularly appropriate when knotting monofilament sutures.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third ( $\frac{1}{3}$ ) to one-half ( $\frac{1}{2}$ ) of the

distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in Sharps containers.

### **ADVERSE REACTIONS**

Adverse effects associated with the use of this device include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal to mild inflammatory tissue reaction, and transitory local irritation at the wound site.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

### **HOW SUPPLIED**

PRONOVA sutures, pigmented, are available as sterile strands in U.S.P. sizes 10-0 through 8-0 (metric sizes 0.2-0.4) and 6-0 through 2 (metric sizes 0.7-5.0).

PRONOVA sutures, clear, are available as sterile strands in U.S.P. sizes 6-0 through 2 (metric sizes 0.7-5.0). Size 7-0 (metric size 0.5) PRONOVA sutures, pigmented and clear are U.S.P. except for diameter. All PRONOVA sutures are available in a variety of lengths, with permanently attached needles.

PRONOVA sutures, pigmented and clear are also available as sterile strands in U.S.P. sizes 0 through 2 (metric sizes 3.5-5.0) attached to CONTROL RELEASE\* removable needles.

PRONOVA sutures are available in one, two, and three dozen boxes.

## ***PROLENE Polypropylene Hernia System***

### **NONABSORBABLE SYNTHETIC SURGICAL MESH**

#### **DESCRIPTION**

The PROLENE polypropylene Hernia System is a sterile, pre-shaped, three-dimensional device constructed of an onlay patch connected by a mesh cylinder to a circular underlay patch. The material is undyed PROLENE polypropylene mesh constructed of knitted nonabsorbable polypropylene filaments.

#### **ACTIONS/PERFORMANCE**

The PROLENE Hernia System is a nonabsorbable mesh used to reinforce or bridge inguinal hernia deficiencies to provide extended support during and following wound healing. Animal studies show that implantation of PROLENE mesh elicits a minimum to slight inflammatory reaction, which is transient and is followed by the deposition of a thin fibrous layer of tissue which can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The mesh remains soft and pliable, and normal wound healing is not noticeably impaired. The material is neither absorbed nor is it subject to degradation or weakening by the action of tissue enzymes.

#### **INDICATIONS**

This product is indicated for the repair of indirect and direct inguinal hernia defects.

#### **WARNINGS**

The PROLENE Hernia System is provided by ETHICON, INC. as a sterile product. This device is for single use only. Do not resterilize. Discard opened packages and unused product. When this device is used in infants or children with future growth potential, the surgeon should be aware that this product will not stretch significantly as the patient grows. The PROLENE Hernia System should only be used in contaminated wounds with the understanding that subsequent infection may require removal of the device.

#### **PRECAUTIONS**

Sutures or clips, if necessary, should be placed such that a minimum of 6.5mm (1/4") of mesh should extend beyond the suture line.

#### **ADVERSE REACTIONS**

Potential adverse reactions are those typically associated with surgically implantable materials which include infection potentiation, inflammation, adhesion formation, fistula formation and extrusion.

#### **STERILITY**

The PROLENE Hernia System is sterilized by Ethylene Oxide. Do not resterilize. Do not use if package is opened or damaged. Discard open, unused product.

#### **STORAGE**

Recommended storage conditions: below 25°C, 77°F, away from moisture and direct heat. Do not use after expire date.

#### **HOW SUPPLIED**

The PROLENE Hernia System is available sterile, undyed in two sizes-medium and large.

#### **CAUTION**

Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

## ***PROLENE Polypropylene Mesh***

### **NONABSORBABLE SYNTHETIC SURGICAL MESH**

#### **DESCRIPTION**

PROLENE polypropylene mesh is constructed of knitted filaments of extruded polypropylene identical in composition to that used in PROLENE Polypropylene Suture, Nonabsorbable Surgical Sutures, U.S.P. (ETHICON, INC.). The mesh is approximately 0.027 inches thick. This material, when used as a suture, has been reported to be nonreactive and to retain its strength indefinitely in clinical use.

PROLENE mesh is knitted by a process which interlinks each fiber junction and which provides for elasticity in both directions. This construction permits the mesh to be cut into any desired shape or size without unraveling. The fiber junctions are not subject to the same work fatigue exhibited by more rigid metallic meshes. This bi-directional elastic property allows adaptation to various stresses encountered in the body.

#### **ACTIONS**

PROLENE mesh is a nonabsorbable mesh used to span and reinforce traumatic or surgical wounds to provide extended support during and following wound healing. Animal studies show that implantation of PROLENE mesh elicits a minimum to slight inflammatory reaction, which is transient and is followed by the deposition of a thin fibrous layer of tissue which can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The mesh remains soft and pliable, and normal wound healing is not noticeably impaired. The material is not absorbed nor is it subject to degradation or weakening by the action of tissue enzymes.

#### **INDICATIONS**

This mesh may be used for the repair of hernia and other fascial deficiencies that require the addition of a reinforcing or bridging material to obtain the desired surgical result.

#### **CONTRAINDICATIONS**

When this mesh is used in infants or children with future growth potential, the surgeon should be aware that this product will not stretch significantly as the patient grows. PROLENE mesh in contaminated wounds should be used with the understanding that subsequent infection may require removal of the material.

#### **WARNINGS**

PROLENE mesh is provided by ETHICON, INC. as a sterile product. Resterilization of the device is NOT recommended. However, testing has demonstrated that reprocessing of unused PROLENE mesh which has been removed from the package will not be adversely affected when exposed not more than one time to conventional steam autoclave conditions at 250°F (121°C) for 20 minutes. Reprocessing under any other condition or by any other means is neither recommended nor endorsed by ETHICON, INC. PROLENE mesh should not be flash autoclaved.

If this product should become stained with blood or soiled, it should not be resterilized for reuse. When reprocessed as outlined above, it is the responsibility of the end-user to assure sterility of the product via a validated sterilization process as ETHICON, INC. has no control over environmental conditions the product may encounter prior to, during, or after reprocessing.

#### **PRECAUTIONS**

A minimum of 6.5mm (1/4") of mesh should extend beyond the suture line.

#### **ADVERSE REACTIONS**

Potential adverse reactions are those typically associated with surgically implantable materials which include infection potentiation, inflammation, adhesion formation, fistula formation and extrusion.

#### **INSTRUCTIONS FOR USE**

It is recommended that nonabsorbable sutures be placed 6.5mm to 12.5 mm (1/4" to 1/2") apart at a distance approximately 6.5mm (1/4") from edge of the mesh. Some surgeons prefer to suture an uncut section of mesh that is considerably larger than the defect into position over the wound. The opposite sides are then sutured to assure proper closure under correct tension. When the margin sutures have all been placed the extra mesh is trimmed away.

#### **HOW SUPPLIED**

PROLENE mesh is available in single packets as sterile, undyed (clear) sheets in seven sizes. The sizes available are 2.5cm x 10cm (1" x 4"), 4.6cm x 10.2cm (1.8" x 4"), 6cm x 11cm (2.5" x 4.5"), 6.1cm x 13.7cm (2.4" x 5.4"), 7.6cm x 15cm (3" x 6"), 15cm x 15cm (6" x 6") and 30cm x 30cm (12" x 12"). Each sheet is approximately 0.7mm (0.027") thick.

## ***VICRYL (Polyglactin 910) Knitted Mesh***

### **DESCRIPTION**

VICRYL (polyglactin 910) knitted mesh is prepared from a synthetic absorbable copolymer of glycolide and lactide, derived respectively from glycolic and lactic acids. This knitted mesh is prepared from uncoated, undyed fiber identical in composition to that used in VICRYL (polyglactin 910) synthetic absorbable suture, which has been found to be inert, nonantigenic, nonpyrogenic and to elicit only a mild tissue reaction during absorption.

VICRYL knitted mesh is intended for use as a buttress to provide temporary support during the healing process.

### **ACTIONS**

Two important characteristics describe the in vivo function and behavior of VICRYL knitted mesh: reinforced wound strength and the rate of absorption (loss of mass).

The dehiscence force of healing abdominal wounds in rats closed with size 4-0 absorbable sutures was compared with corresponding wounds closed with size 4-0 absorbable sutures and reinforced with VICRYL knitted mesh. In this animal model, the strength of the incision, when supported by the mesh, was significantly greater than the sutured incisional wound. Explanted VICRYL knitted mesh, which, before implantation had an initial average burst strength of 63 lbs., was found to have 80% of its original burst strength remaining after fourteen days in vivo.

Subcutaneous implantation studies in rats indicate that the absorption of VICRYL mesh material is minimal until about six weeks post implantation and essentially complete between 60 and 90 days.

### **INDICATIONS**

VICRYL knitted mesh may be used wherever temporary wound or organ support is required, particularly in instances in which compliant and stretchable support material is desired and containment of wound transudate is not required. VICRYL knitted mesh may be cut to the shape or size desired for each specific application.

### **CONTRAINDICATIONS**

Because VICRYL knitted mesh is absorbable, it should not be used where extended wound or organ support is required.

### **WARNINGS**

**DO NOT RESTERILIZE.**

The safety and effectiveness of VICRYL knitted mesh in neural and in cardiovascular tissue has not been established.

### **PRECAUTIONS**

None.

### **ADVERSE REACTIONS**

No significant clinical adverse reactions to the mesh have been reported.

### **DIRECTIONS FOR USE**

It is recommended that absorbable or nonabsorbable sutures be placed  $\frac{1}{4}$  to  $\frac{1}{2}$  inch (6 to 12mm) apart at a distance approximately  $\frac{1}{4}$  inch (6mm) from the edge of the mesh. Some surgeons prefer to suture a mesh larger than the defect into position over the defect. The edges are then sutured to assure proper closure under correct tension. When all margin sutures have been placed, the excess mesh is trimmed away, leaving at least  $\frac{1}{4}$  inch of mesh extending beyond the suture line.

### **HOW SUPPLIED**

VICRYL knitted mesh is available in single packets as a sterile, undyed fabric mesh in single sheet sizes of approximately 6 x 6 inches and 12 x 12 inches (15 x 15 centimeters and 30 x 30 centimeters).

## ***VICRYL (Polyglactin 910) Woven Mesh***

### **DESCRIPTION**

VICRYL (polyglactin 910) woven mesh is prepared from a synthetic absorbable copolymer of glycolide and lactide, derived respectively from glycolic and lactic acids. This tightly woven mesh is prepared from uncoated, undyed fiber identical in composition to that used in VICRYL (polyglactin 910) synthetic absorbable suture, which has been found to be inert, nonantigenic, nonpyrogenic and to elicit only a mild tissue reaction during absorption.

VICRYL woven mesh is intended for use as a buttress to provide temporary support during the healing process.

### **ACTIONS**

Two important characteristics describe the in vivo function and behavior of VICRYL woven mesh: reinforced wound strength and the rate of absorption (loss of mass).

The dehiscence force of healing abdominal wounds in rats closed with size 4-0 absorbable sutures was compared with corresponding wounds closed with size 4-0 absorbable sutures and reinforced with VICRYL woven mesh. In this animal model, the strength of the incision, when supported by the mesh, was significantly greater than the sutured incisional wound. Explanted VICRYL woven mesh, which, before implantation had an initial average burst strength of approximately 121 lbs., was found to have approximately 23% of its original burst strength remaining after fourteen days in vivo.

Subcutaneous implantation studies in rats indicate that the absorption of VICRYL mesh material is minimal until about six weeks post implantation and essentially complete between 60 and 90 days.

### **INDICATIONS**

VICRYL woven mesh may be used wherever temporary wound or organ support is required. The woven mesh structure is less porous than VICRYL knitted mesh. It is indicated in instances in which containment of wound transudate is desirable. VICRYL woven mesh may be cut to the shape or size desired for each specific application.

### **CONTRAINDICATIONS**

Because VICRYL woven mesh is absorbable, it should not be used where extended wound or organ support is required.

### **WARNINGS**

#### **DO NOT RESTERILIZE.**

The safety and effectiveness of VICRYL woven mesh in neural tissue and in cardiovascular tissue has not been established.

### **PRECAUTIONS**

None.

### **ADVERSE REACTIONS**

None known.

### **DIRECTIONS FOR USE**

It is recommended that absorbable or nonabsorbable sutures be placed  $\frac{1}{4}$  to  $\frac{1}{2}$  inch (6 to 12mm) apart at a distance approximately  $\frac{1}{4}$  inch (6mm) from the edge of the mesh. Some surgeons prefer to suture a mesh larger than the defect into position over the defect. The edges are then sutured to assure proper closure under correct tension. When all margin sutures have been placed, the excess mesh is trimmed away, leaving at least  $\frac{1}{4}$  inch of mesh extending beyond the suture line.

### **HOW SUPPLIED**

VICRYL woven mesh is available in single packets as a sterile, undyed fabric mesh in single sheet sizes of approximately 6 x 6 inches and 12 x 12 inches (15 x 15 centimeters and 30 x 30 centimeters).

## ***MERSILENE Polyester Fiber Mesh***

NONABSORBABLE SYNTHETIC

SURGICAL MESH, STERILE

### **DESCRIPTION**

MERSILENE Polyester Fiber Mesh is constructed from polyethylene terephthalate, the same material used to make MERSILENE Polyester Fiber Suture, Nonabsorbable Surgical Suture, U.S.P. (ETHICON, INC.). MERSILENE Polyester Fiber Mesh affords the excellent strength, durability and surgical adaptability, along with maximal porosity for necessary tissue ingrowth. The mesh is approximately 0.010 inches thick and is a highly flexible and compliant material.

MERSILENE mesh is knitted by a process which interlinks each fiber junction and which provides for elasticity in both directions. This construction permits the mesh to be cut into any desired shape or size without unraveling. The fiber junctions are not subject to the same work fatigue exhibited by more rigid metallic meshes. This bi-directional elastic property allows adaption to various stresses encountered in the body.

### **ACTIONS**

MERSILENE mesh is a nonabsorbable mesh used to span and reinforce traumatic or surgical wounds to provide extended support during and following wound healing. Animal studies show that implantation of MERSILENE mesh elicits a minimum to slight inflammatory reactions, which is transient and is followed by the deposition of a thin fibrous layer of tissue which can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The mesh remains soft and pliable, and normal wound healing is not noticeably impaired. The material is not absorbed nor is it subject to degradation or weakening by the action of tissue enzymes.

### **INDICATIONS**

This mesh may be used for the repair of hernia and other fascial deficiencies that require the addition of a reinforcing or bridging material to obtain the desired surgical result.

### **CONTRAINDICATIONS**

When this mesh is used in infants or children with future growth potential, the surgeon should be aware that this product will not stretch significantly as the patient grows.

MERSILENE polyester fiber mesh in contaminated wounds should be used with the understanding that subsequent infection may require removal of the material.

### **WARNINGS**

MERSILENE mesh is provided by ETHICON, INC. as a sterile product. Unused MERSILENE mesh which has been removed from the package may be resterilized not more than one time by a conventional steam autoclaving process at conditions of 250°F (121°C) for 20 minutes. MERSILENE mesh may also be flash autoclaved not more than one time at conditions of 270°F (132°C) for 10 minutes. Resterilization under any other conditions or by any other means is neither recommended nor endorsed by ETHICON, INC.

If this product should become stained with blood or soiled, it should not be resterilized for reuse.

### **PRECAUTIONS**

A minimum of 6.5mm (1/4 inch) of mesh should extend beyond the suture line.

### **ADVERSE REACTIONS**

No significant adverse clinical reactions to MERSILENE mesh have been reported. The use of nonabsorbable MERSILENE mesh in a wound that is contaminated or infected could lead to fistula formation and/or extrusion of the mesh.

### **INDICATIONS FOR USE**

It is recommended that nonabsorbable sutures be placed 6.5 to 12.5mm (1/4 to 1/2 inch) apart at a distance approximately 6.5mm (1/4 inch) from edge of the mesh. Some surgeons prefer to suture an uncut section of mesh that is considerably larger than the defect into position over the wound. The opposite sides are then sutured to assure proper closure under correct tension. When the margin sutures have all been placed, the extra mesh is trimmed away.

### **HOW SUPPLIED**

MERSILENE mesh is available in single packets as sterile, undyed (white) sheets in two sizes. The sizes available are 6 x 11cm (2.5 x 4.5 inches) and 30 x 30cm (12 x 12 inches). Each sheet is 0.25mm (0.010 inch) thick.

## ***DERMABOND Topical Skin Adhesive page 1 of 3***

(2 - OCTYL CYANOACRYLATE)

### **DESCRIPTION**

DERMABOND Topical Skin Adhesive is a sterile, liquid topical skin adhesive containing a monomeric (2-octyl cyanoacrylate) formulation and the colorant D & C Violet #2. It is provided in a single use applicator packaged in a blister pouch. The applicator is comprised of a crushable glass ampule contained within a plastic vial with attached applicator tip. As applied to the skin, the liquid adhesive is slightly more viscous than water and polymerizes within minutes. See DIRECTIONS FOR USE.

### **INDICATIONS**

DERMABOND adhesive is intended for topical application only to hold closed easily approximated skin edges of wounds from surgical incisions, including punctures from minimally invasive surgery, and simple, thoroughly cleansed, trauma-induced lacerations. DERMABOND adhesive may be used in conjunction with, but not in place of, subcuticular sutures.

### **CONTRAINDICATIONS**

- DERMABOND adhesive is contraindicated for use on any wounds with evidence of active infection, gangrene, or wounds from decubitus etiology.
- DERMABOND adhesive is contraindicated for use on mucosal surfaces or across mucocutaneous junctions (e.g., oral cavity, lips), or on skin which may be regularly exposed to body fluids or with dense natural hair, (e.g., scalp).
- DERMABOND adhesive is contraindicated for use on patients with a known hypersensitivity to cyanoacrylate or formaldehyde.

### **WARNINGS**

- DERMABOND adhesive is a fast setting adhesive capable of adhering to most body tissue and many other materials, such as latex gloves and stainless steel. Inadvertent contact with any body tissue, and any surfaces or equipment that are not disposable or that cannot be readily cleaned with a solvent such as acetone should be avoided.
- Polymerization of DERMABOND adhesive may be accelerated by water or fluids containing alcohol: DERMABOND adhesive should not be applied to wet wounds.
- DERMABOND adhesive should not be applied to the eye. If contact with the eye occurs, flush the eye copiously with saline or water. If residual adhesive remains, apply topical ophthalmic ointment to help loosen the bond and contact an ophthalmologist.
- DERMABOND adhesive should not be used below the skin because the polymerized material is not absorbed by tissue and can elicit a foreign body reaction.
- DERMABOND adhesive should not be used in high skin tension areas or across areas of increased skin tension, such as knuckles, elbows, or knees, unless the joint will be immobilized during the skin healing period.

- DERMABOND adhesive treated wounds should be monitored for signs of infection. Wounds with signs of infection, such as erythema, edema, warmth, pain and pus, should be evaluated and treated according to standard practice for infection.
- DERMABOND adhesive should not be used on wound sites that will be subjected to repeated or prolonged moisture or friction.
- DERMABOND adhesive should only be used after wounds have been cleaned and debrided in accordance with standard surgical practice. Local anesthetic should be used when necessary to assure adequate cleansing and debridement.
- Excessive pressure of the applicator tip against the wound edges or surrounding skin can force the wound edges apart and allow adhesive into the wound. Adhesive within the wound could delay wound healing and/or result in adverse cosmetic outcome. Therefore, DERMABOND adhesive should be applied with a very light brushing motion of the applicator tip over easily approximated wound edges.
- DERMABOND adhesive polymerizes through an exothermic reaction in which a small amount of heat is released. With the proper technique of applying DERMABOND adhesive in multiple thin layers (at least three) onto a dry wound and allowing time for polymerization between applications, heat is released slowly and the sensation of heat or pain experienced by the patient is minimized. However, if DERMABOND adhesive is applied so that large droplets of liquid are allowed to remain unspread, the patient may experience a sensation of heat or discomfort.
- DERMABOND adhesive is packaged for single patient use. Discard remaining opened material after each wound closure procedure.
- Do not resterilize DERMABOND adhesive.
- Do not place DERMABOND adhesive in a procedure pack/tray that is to be sterilized prior to use. Exposure of DERMABOND adhesive, after its final manufacture, to excessive heat (as in autoclaves or ethylene oxide sterilization) or radiation (such as gamma or electron beam), is known to increase its viscosity and may render the product unusable.

#### **PRECAUTIONS**

- Do not apply liquid or ointment medications or other substances to the wound after closure with DERMABOND adhesive, as these substances can weaken the polymerized film and allow for wound dehiscence. DERMABOND adhesive permeability by topical medications has not been studied.
- DERMABOND adhesive permeability by fluids is not known and has not been studied.

- DERMABOND adhesive is a free flowing liquid slightly more viscous than water. To prevent inadvertent flow of liquid DERMABOND adhesive to unintended areas: (1) the wound should be held in a horizontal position, with DERMABOND adhesive applied from above, and (2) DERMABOND adhesive should be applied in multiple (at least 3), thin layers rather than in a few large droplets.
- DERMABOND adhesive should be used immediately after crushing the glass ampule as the liquid adhesive will not flow freely from the applicator tip after a few minutes.
- If unintended bonding of intact skin occurs, peel, but do not pull the skin apart. Petroleum jelly or acetone may help loosen the bond. Other agents such as water, saline, betadine, Hibiclens, or soap, are not expected to immediately loosen the bond.
- Safety and effectiveness of DERMABOND adhesive on wounds of patients with peripheral vascular disease, insulin dependent diabetes mellitus, blood clotting disorders, personal or family history of keloid formation or hypertrophy, burst stellate lacerations, have not been studied.
- Safety and effectiveness of DERMABOND adhesive on the following wounds have not been studied: animal or human bites, puncture or stab wounds.
- Safety and effectiveness on wounds that have been treated with DERMABOND adhesive and then exposed for prolonged periods to direct sunlight or tanning lamps have not been studied.
- Safety and effectiveness of DERMABOND adhesive on wounds in vermilion surfaces has not been studied.

## ADVERSE REACTIONS

- Adverse reactions encountered during clinical study:

<b>Clinical Study Outcomes</b>	<b>No Subcuticular Sutures</b>		<b>With Subcuticular Sutures</b>	
	DERMABOND	Control	DERMABOND	CONTROL
	N(%)	N(%)	N(%)	N(%)
<b>Accounting</b>				
N, patients enrolled	240	243	167	168
N, patients treated	239	242	167	166
Patients completed	228 (95%)	215 (88%)	164 (98%)	162 (96%)
<b>Adverse Reactions</b>				
Suspected infection*	8 (3.6%)	2 (0.9%)	6 (3.6%)	2 (1.2%)
Wound type #Lacerations	8	2	1	0
#Incisions	0	0	5	2
<b>Dehiscence with Need for Retreatment</b>	6 (2.5%)	5 (2.1%)	3 (1.8%)	0
<b>Accute Inflammation</b>				
Erythema				

## ***DERMABOND Topical Skin Adhesive page 3 of 3***

Follow-up was at 5-10 days and at 3 months. All wounds were assessed by visual inspection at 5-10 days after wound closure. The types of wounds treated in the study were 46.1% lacerations and 53.9% incisions. The incisions were comprised of 47.8% excisions of skin lesions, 27.3% minimally invasive surgery punctures, and 24.8% general surgery incisions.

For wounds closed without subcuticular stitches, mean wound length was 1.5cm, mean wound width was 2.5mm, and mean wound depth was 5.8mm. For wounds closed with subcuticular stitches, mean wound length was 3.2cm, mean wound width was 5.3mm, and mean wound depth was 3.8mm.

If the primary method of closure was insufficient for closure, an additional securing device was placed. The time to perform treatment included the time required later to remove the closure device when applicable.

The Modified Hollander Cosmesis Scale (MHCS), a validated scale, was used to evaluate cosmesis at three months: step-off borders, edge inversion, contour irregularities, excess inflammation, wound margin separation, and overall appearance.

### **DIRECTIONS FOR USE**

1. The application of DERMABOND adhesive requires thorough wound cleansing. Follow standard surgical practice for wound preparation before application of DERMABOND adhesive (i.e., anesthetize, irrigate, debride, obtain hemostasis and close deep layers).
2. Pat the wound dry with dry, sterile gauze to assure direct tissue contact for adherence of DERMABOND adhesive to the skin. Moisture accelerates the polymerization of DERMABOND adhesive and may affect wound closure results.
3. To prevent inadvertent flow of liquid DERMABOND adhesive to unintended areas of the body, the wound should be held in a horizontal position and DERMABOND adhesive should be applied from above the wound.
4. DERMABOND adhesive should be used immediately after crushing the glass ampule, since the liquid adhesive will flow freely from the tip for only a few minutes. Remove the applicator from the blister pouch. Hold the applicator with the thumb and a finger and away from the patient to prevent any unintentional placement of the liquid DERMABOND adhesive into the wound or on the patient. While holding the applicator, and with applicator tip pointed upward, apply pressure at the midpoint of the ampule to crush the inner glass ampule. Invert and gently squeeze the applicator just sufficiently to express the liquid DERMABOND adhesive to moisten the applicator tip.
5. Approximate wound edges with gloved fingers or sterile forceps. Slowly apply the liquid DERMABOND adhesive in multiple (at least 3) thin layers to the surface of the approximated wound edges using a gentle brushing motion. Wait approximately 30 seconds between applications or layers. Maintain manual approximation of the wound edges for approximately 60 seconds after the final layer.

**NOTE:** DERMABOND adhesive polymerizes through an exothermic reaction. If the liquid DERMABOND adhesive is applied so that large droplets are allowed to remain without being evenly spread, the patient may experience a sensation of heat or discomfort. The sensation may be higher on sensitive tissues. This can be minimized by applying DERMABOND adhesive in multiple thin layers (at least 3).

**NOTE:** Excessive pressure of the applicator tip against the wound edges or surrounding skin can result in forcing the wound edges apart and allowing DERMABOND adhesive into the wound. DERMABOND adhesive within the wound could delay wound healing and/or result in adverse cosmetic outcome.

**NOTE:** Full apposition strength is expected to be achieved about 2.5 minutes after the final layer is applied, although the top adhesive layer may remain tacky for up to approximately 5 minutes. Full polymerization is expected when the top DERMABOND adhesive layer is no longer sticky.

6. Do not apply liquid or ointment medications onto wounds closed with DERMABOND adhesive because these substances can weaken the polymerized film, leading to wound dehiscence.

7. Protective dry dressings, such as gauze, may be applied only after DERMABOND adhesive film is completely solid/polymerized: not tacky to the touch (approximately five minutes after application). Allow the top layer to fully polymerize before applying a bandage.

If a dressing, bandage, adhesive backing or tape is applied before complete polymerization, the dressing can adhere to the film. When the dressing is removed, the film can be disrupted from the skin and wound dehiscence can occur.

8. Patients should be instructed to not pick at the polymerized film of DERMABOND adhesive. Picking at the film can disrupt its adhesion to the skin and cause dehiscence of the wound. Picking at the film can be discouraged by an overlying dressing.

9. Apply a dry protective dressing for children or other patients who may not be able to follow instructions for proper wound care.

10. Patients treated with DERMABOND adhesive should be provided the printed instruction sheet entitled How to Care For Your Wound After It's Treated With DERMABOND Topical Skin Adhesive. This instruction sheet should be reviewed with each patient or guardian to assure understanding of the proper care for the treatment site.

11. Patients should be instructed that until the polymerized film of DERMABOND adhesive has sloughed naturally (usually in 5-10 days), there should be only transient wetting of the treatment site. Patients may shower and bathe the site gently. The site should not be scrubbed, soaked, or exposed to prolonged wetness until after the film has sloughed naturally and the wound has healed closed. Patients should be instructed not to go swimming during this period.

12. If removal of DERMABOND adhesive is necessary for any reason, carefully apply petroleum jelly or acetone to the DERMABOND film to help loosen the bond. Peel off the film, do not pull the skin apart.

#### **HOW SUPPLIED**

DERMABOND adhesive is supplied sterile, in a pre-filled, single-use applicator. The applicator is comprised of a crushable glass ampule contained with a plastic vial with attached applicator tip. The applicator contains 0.5 ml of liquid adhesive. The applicator is packaged in a blister pouch to maintain the device sterile until opened or damaged.

DERMABOND adhesive is available in boxes of 12 applicators.

#### **STORAGE**

Recommended storage conditions: below 30°C, 86°F, away from moisture and direct heat. Do not use after expiry date.

#### **STERILITY**

DERMABOND adhesive is originally sterilized by dry heat and ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged. Discard any unused material following completion of medical procedure.

#### **STERILE SINGLE USE ONLY**

#### **REPORTING**

Physicians should use the following toll free number 1-800-255-2500 (valid in U.S.A. only), when reporting adverse reactions or potentially threatening complications involving DERMABOND adhesive.

#### **CAUTION**

Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Manufactured for ETHICON, INC. by Closure Medical Corp.

## ***SURGICEL, SURGICEL Fibrillar, and SURGICEL Nu-Knit Absorbable Hemostats page 1 of 2***

(OXIDIZED REGENERATED CELLULOSE)

FOR SURGICAL USE

### **BRIEF SUMMARY**

See package insert for complete prescribing information and endoscopic application.

### **DESCRIPTION**

SURGICEL absorbable hemostat is a sterile absorbable knitted fabric prepared by the controlled oxidation of regenerated cellulose. The fabric is white with a pale yellow cast and has a faint, caramel-like aroma. It is strong and can be sutured or cut without fraying. It is stable and should be stored at controlled room temperature. A slight discoloration may occur with age, but this does not affect performance.

The fibrillar form of the product allows the surgeon to grasp with forceps any amount of SURGICEL fibrillar hemostat needed to achieve hemostasis at a particular bleeding site. The fibrillar form may be more convenient than the knitted form for hard to reach or irregularly shaped bleeding sites. Although it is easy to pull the desired amount of SURGICEL fibrillar hemostat from the entire supply, the group of selected fibers continue to cohere to one another and application to the bleeding site is easily controlled. Unwanted dispersal over the operative site does not occur.

### **INDICATIONS**

SURGICEL hemostat (oxidized regenerated cellulose) is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial hemorrhage when ligation or other conventional methods of control are impractical or ineffective. SURGICEL, SURGICEL fibrillar, and SURGICEL NU-KNIT hemostats can be cut to size for use in endoscopic procedures.

### **CONTRAINDICATIONS**

Although packing or wadding sometimes is medically necessary, SURGICEL hemostat should not be used in this manner, unless it is to be removed after hemostasis is achieved (See WARNINGS and PRECAUTIONS).

SURGICEL hemostat should not be used for implantation in bone defects, such as fractures, since there is a possibility of interference with callus formation and a theoretical chance of cyst formation.

When SURGICEL hemostat is used to help achieve hemostasis in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, or the optic nerve and chiasm, it must always be removed after hemostasis is achieved since it will swell and could exert unwanted pressure.

SURGICEL hemostat should not be used to control hemorrhage from large arteries.

SURGICEL hemostat should not be used on non-hemorrhagic serous oozing surfaces, since body fluids other than whole blood, such as serum, do not react with SURGICEL hemostat to produce satisfactory hemostatic effect.

SURGICEL hemostat is an absorbable hemostat, and should not be used as an adhesion prevention product.

### **WARNINGS**

SURGICEL hemostat is supplied sterile and as the material is not compatible with autoclaving or ethylene oxide sterilization, SURGICEL hemostat should not be resterilized.

SURGICEL hemostat is not intended as a substitute for careful surgery and the proper use of sutures and ligatures.

Closing SURGICEL hemostat in a contaminated wound without drainage may lead to complications and should be avoided.

The hemostatic effect of SURGICEL hemostat is greater when it is applied dry; therefore it should not be moistened with water or saline.

SURGICEL hemostat should not be impregnated with anti-infective agents or with any other materials such as buffering or hemostatic substances. Its hemostatic effect is not enhanced by the addition of thrombin, the activity of which is destroyed by the low pH of the product.

Although SURGICEL hemostat may be left in situ when necessary, it is advisable to remove it once hemostasis is achieved. It must always be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm regardless of the type of surgical procedure because SURGICEL hemostat, by swelling, may exert pressure resulting in paralysis and/or nerve damage. Dislodgement of SURGICEL hemostat could possible occur by means such as repacking, further intraoperative manipulation, lavage, exaggerated respiration, etc. There

have been reports that in procedures such as lobectomy, laminectomy and repair of frontal skull fracture and lacerated lobe that SURGICEL hemostat, when left in the patient after closure, migrated from the site of application into foramina in bone around the spinal cord resulting in paralysis and, in another case, the left orbit of the eye, causing blindness. While these reports cannot be confirmed, special care must be taken by physicians, *regardless of the type of surgical procedure*, to consider the advisability of removing SURGICEL hemostat after hemostasis is achieved.

Although SURGICEL hemostat is bactericidal against a wide range of pathogenic microorganisms, it is not intended as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or prevent postoperative infections.

#### **PRECAUTIONS**

Use only as much SURGICEL hemostat as is necessary for hemostasis, holding it firmly in place until bleeding stops. Remove any excess before surgical closure in order to facilitate absorption and minimize the possibility of foreign body reaction.

In urological procedures, minimal amounts of SURGICEL hemostat should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product. Since absorption of SURGICEL hemostat could be prevented in chemically cauterized areas, its use should not be preceded by application of silver nitrate or any other escharotic chemicals.

## ***SURGICEL, SURGICEL Fibrillar, and SURGICEL Nu-Knit Absorbable Hemostats page 2 of 2***

If SURGICEL hemostat is used temporarily to line the cavity of large open wounds, it should be placed so as not to overlap the skin edges. It should also be removed from open wounds by forceps or by irrigation with sterile water or saline solution after bleeding has stopped.

Precautions should be taken in otorhinolaryngologic surgery to assure that none of the material is aspirated by the patient. (Examples: controlling hemorrhage after tonsillectomy and controlling epistaxis.)

Care should be taken not to SURGICEL hemostat too tightly when it is used as a wrap during vascular surgery (see ADVERSE REACTIONS).

Endoscopic procedures should be performed only by persons having adequate training and familiarity with endoscopic techniques. Consult medical literature relative to techniques, complications and hazards prior to performance of any endoscopic procedure.

A thorough understanding of the principles and techniques involved in laparoscopic laser and electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Refer to appropriate electrosurgical system users manual for use indications and instructions to ensure that all safety precautions are followed. When endoscopic instruments and accessories from different manufacturers are employed together during a procedure, verify their compatibility prior to initiation of the procedure and ensure that isolation or grounding is not compromised.

### **ADVERSE REACTIONS**

"Encapsulation" of fluid and foreign body reactions have been reported.

There have been reports of stenotic effect when SURGICEL hemostat has been applied as a wrap during vascular surgery. Although it has not been established that the stenosis was directly related to the use of SURGICEL hemostat, it is important to be cautious and avoid applying the material tightly as a wrapping.

Paralysis and nerve damage have been reported when SURGICEL hemostat was used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm.

While most of these reports have been in connection with laminectomy, reports of paralysis have also been received in connection with other procedures. Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when SURGICEL hemostat was placed in the anterior cranial fossa (See WARNINGS and PRECAUTIONS).

Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy have been reported. There has been one report of a blocked ureter after kidney resection, in which postoperative catheterization was required.

Occasional reports of "burning" and "stinging" sensations and sneezing when SURGICEL hemostat has been used as packing in epistaxis, are believed due to the low pH of the product.

Burning has been reported when SURGICEL hemostat was applied after nasal polyp removal and after hemorrhoidectomy. Headache, burning, stinging, and sneezing in epistaxis and other rhinological procedures, and stinging when SURGICEL hemostat was applied on surface wounds (varicose ulcerations, dermabrasions, and donor sites) also have been reported.

### **HOW SUPPLIED**

Sterile SURGICEL fibrillar absorbable hemostat is supplied as stable fiber in an envelopes in the following sizes.

Code No. 1961 1 in. x 2 in. (2.5cm x 5.1cm)

Code No. 1962 2 in. x 4 in. (5.1cm x 10.2cm)

Code No. 1963 4 in x 4 in. (10.2cm x 10.2cm).

Sterile SURGICEL absorbable hemostat (oxidized regenerated cellulose) is supplied as knitted fabric strips in envelopes in the following sizes.

Code No. 1951 2 in. x 14 in. (5.1cm x 35.6cm)

Code No. 1952 4 in. x 8 in. (10.2cm x 20.3cm)

Code No. 1953 2 in. x 3 in. (5.1cm x 7.6cm)

Code No. 1955 1<sup>1</sup>/<sub>2</sub> in. x 2 in. (1.3cm x 5.1cm).

SURGICEL NU-KNIT absorbable hemostat

Code No. 1940 1 in. x 1 in. (2.5cm x 2.5cm)

Code No. 1941 1 in. x 3.5 in. (2.5cm x 8.9cm)

Code No. 1943 3 in. x 4 in. (7.6cm x 10.2cm)

Code No. 1946 6 in. x 9 in. (15.2cm x 22.9cm).

**STORAGE**

Store at controlled temperature 15?-30?C (59?-86?F).

**CAUTION**

Federal law restricts this device to sale by or on the order of a physician.

**DESCRIPTION**

INSTAT collagen absorbable hemostat is a purified and lyophilized bovine dermal collagen. The material, prepared as a sponge-like pad, is lightly cross-linked, sterile, non-pyrogenic, and absorbable. Hemostatic activity, which is an inherent property of collagen, is largely dependent on the basic helical structure of this protein. The helical structure of native collagen is preserved during the manufacture of INSTAT hemostat. When collagen comes into contact with blood, platelets aggregate on the collagen and release coagulation factors which, together with plasma factors result in the formation of fibrin, and finally in the formation of a clot.

**INDICATIONS**

INSTAT hemostat is indicated in surgical procedures (other than in neurosurgical and ophthalmological surgery) for use as an adjunct to hemostasis when control of bleeding by ligature or other conventional methods is ineffective or impractical.

**CONTRAINDICATIONS**

INSTAT hemostat should not be used in the closure of skin incisions as it may interfere with the healing of skin edges. This interference is due to simple mechanical interposition of dry collagen and not due to any intrinsic interference with wound healing. It has been reported with another absorbable collagen hemostat that, in filling porosities of cancellous bone, collagen may reduce the bonding strength of methylmethacrylate. Therefore, INSTAT hemostat should not be applied on bone surfaces to which prosthetic materials are to be attached with methylmethacrylate adhesives.

**WARNINGS**

INSTAT hemostat is inactivated by autoclaving. It should not be resterilized. As with any foreign substance, use in contaminated wounds may enhance infection.

INSTAT hemostat should not be used in instances of pumping arterial hemorrhage.

INSTAT hemostat should not be used where blood or other fluids have pooled or in cases where the point of hemorrhage is submerged. INSTAT hemostat will not act as a tampon or plug in a bleeding site nor will it close off an area of blood collection behind a tampon.

Only the amount of INSTAT hemostat necessary to provide hemostasis should be used. The long-term effects of leaving INSTAT hemostat in situ are unknown. Opened, unused INSTAT hemostat should be discarded because it cannot be resterilized.

**PRECAUTIONS**

As with other hemostatic agents, it is not recommended that INSTAT hemostat be left in an infected or contaminated space, nor is it recommended for use in persons known to be sensitive to materials of bovine origin. When placed into cavities or closed spaces, care should be exercised to avoid overpacking INSTAT hemostat as it may absorb fluid and expand and press against neighboring structures. In urological procedures, INSTAT hemostat should not be left in the renal pelvis or ureters to eliminate the potential foci for calculus formation.

Safety of this product has not been established in children and pregnant women; therefore, INSTAT hemostat should only be used when benefit to risk clearly warrants its use.

INSTAT hemostat is not intended to be used to treat systemic coagulation disorders.

**ADVERSE REACTIONS**

INSTAT hemostat is a collagen product. Although several types of post-operative complications were observed in INSTAT hemostat treated patients, none were attributed to INSTAT hemostat except one case of fibrotic reaction where INSTAT hemostat involvement could not be ruled out. Adverse reactions reported for other collagen hemostats include hematoma, potentiation of infection, wound dehiscence, inflammation and edema. Other reported adverse reactions that may be related to the use of collagen hemostats include adhesion formation, allergic reaction, foreign body reaction and subgaleal seroma (in a single case). The use of microfibrillar collagen in dental extraction sockets has been reported to increase the incidence of alveolgia. The possibility that all of the above reactions may occur with INSTAT hemostat cannot be excluded.

**ADMINISTRATION**

INSTAT hemostat is applied directly to the bleeding surface with pressure. INSTAT hemostat can be cut to size. The amount needed and the period of time necessary to apply pressure will vary with the type and amount of bleeding to be controlled. Hemostasis time depends upon the type of surgery and degree of pretreatment bleeding. It usually occurred between 2 to 5 minutes with INSTAT hemostat.

INSTAT hemostat maintains its integrity in the presence of blood and is not dispersed when wet. It is easily removed from the site following hemostasis. It is most effective when used dry.

INSTANT hemostat may be left in situ whenever necessary. However, the surgeon, at his discretion, should remove any excess of INSTANT hemostat prior to wound closure. Animal implant studies have demonstrated that absorption and tissue reaction to INSTANT hemostat are similar to those observed with another absorbable collagen hemostatic agent. In these studies, on visual examination, most of INSTANT hemostat was found to be absorbed in 8 to 10 weeks after implantation.

#### **CLINICAL STUDIES**

The safety, effectiveness and handling characteristics of INSTANT hemostat were evaluated in a variety of surgical procedures. The median time to hemostasis for INSTANT hemostat was 3 minutes. Passive Hemagglutination Assay (PHA) and Enzyme-Linked Immunoabsorbent Assay (ELISA) methods have been used to evaluate the immunologic potential for INSTANT hemostat to produce antibodies in patients. These assays revealed mild elevation of antibody titers in both INSTANT hemostat treated patients and patients treated with a collagen control hemostat, confirming that INSTANT hemostat, like other collagen hemostats, is a weak antigen.

#### **HOW SUPPLIED**

INSTANT hemostat is supplied in a sponge-like form in peelable plastic envelopes in the following sizes:

Code No. 1981 1 in. x 2 in. (2.5cm x 5.1cm)

Code No. 1983 3 in. x 4 in. (7.6cm x 10.2cm).

The sterility of the product is guaranteed unless the individual envelope is damaged or opened.

#### **STORAGE**

Store at controlled room temperature 59°-86°F (15°-30°C).

#### **CAUTION**

Federal law restricts this device to sale, distribution, and use by or on the order of a physician.

## **INSTAT MCH Microfibrillar Collagen Hemostat**

### **DESCRIPTION**

INSTAT MCH Microfibrillar Collagen Hemostat is a dry, white absorbent hemostatic agent in a microfibrillar form. The raw material for INSTAT MCH is obtained from bovine deep flexor tendon (achilles tendon), known to be one of the purest sources of collagen that can be readily processed in commercial amounts. Further purification during processing ensures that INSTAT MCH is a consistent material with uniform behavior.

The microfibrillar form of the product allows the surgeon to grasp with forceps any amount of INSTAT MCH needed to achieve hemostasis at a particular bleeding site. The microfibrillar form may be more convenient than the sponge form for hard to reach or irregularly shaped bleeding sites. Although it is easy to pull the desired amount of INSTAT MCH from the entire supply, the group of selected fibers continue to cohere to one another and application to the bleeding site is easily controlled. Unwanted dispersal over the operative site does not occur.

### **INDICATIONS**

INSTAT MCH is recommended for use in surgical procedures (except in Neurological, Urological, and Ophthalmological surgery) as an adjunct to hemostasis when control of bleeding by ligature or conventional procedures is ineffective or impractical.

### **INFORMATION FOR USE**

In contact with blood, collagen is known to cause aggregation of platelets. Platelets deposit in large numbers on the collagen structure, disintegrate and release coagulation factors that, together with plasma factors, enable the formation of fibrin. The microfibrillar structure of INSTAT MCH provides for the additional strengthening of the clot.

Topical application of INSTAT MCH effectively controls bleeding, usually within 2 to 4 minutes, when applied directly to the bleeding site. INSTAT MCH is designed to be completely absorbable if left in situ after hemostasis. In contact with blood, the fibers expand to become a coherent gelatinous mass that conforms to the shape of the bleeding area. If desired, recovery of this mass is easily accomplished.

Absorption of a collagen hemostatic agent was evaluated after subcutaneous and intrahepatic implantation in rats. In 1 out of 5 animals, complete subcutaneous absorption was observed by day 14; and by day 56, 3 out of 4 had complete absorption. Complete intraperitoneal absorption was not observed by day 56. Tissue reactions elicited by the collagen hemostatic agent during implantation are similar to tissue reactions caused by other hemostatic agents.

The collagen hemostatic agent has been evaluated *in-vitro* for the enhancement of bacterial growth of *Staphylococcus aureus* and *Escherichia coli*. Enhancement of bacterial growth did not occur for either organism.

*In vivo* studies using guinea pigs showed that incidence of infection (abscess) of incision sites inoculated with *Staphylococcus aureus* was not enhanced by the presence of the collagen material when compared to another collagen hemostatic agent. However, extent of wound infection tended to be greater than control with the collagen material as was the case with another collagen hemostatic agent tested. This tendency is observed with many foreign substances.

The collagen hemostatic agent was evaluated for potential allergenic sensitivity. A guinea pig maximization study showed that the collagen material did not produce irritation or contact sensitization. A chemical assay of the collagen material compared to one other collagen hemostat showed significantly less specific glycoprotein immunoreactive substances in the collagen material. A hemoagglutination study was conducted evaluating the collagen material as the antigen. There was no agglutination observed.

### **PRECAUTIONS**

As with other hemostatic agents, it is not recommended that INSTAT MCH be left in an infected or contaminated space. Safety of this product has not been established in pregnant women; therefore, it should only be used when benefit to risk clearly warrants its use.

Another microfibrillar collagen agent has been reported to cause interference with the healing of skin edges when used in the closure of skin incisions, and to reduce the strength of methyl methacrylate adhesives used to attach prosthetic devices to bone surfaces.

INSTAT MCH is not intended to be used to treat systemic coagulation disorders.

The long term effects of leaving INSTAT MCH *in situ* are unknown.

### **ADMINISTRATION**

Using forceps, the desired amount of INSTAT MCH can be easily selected from the entire supply. The coherent fibers should then be applied directly to the bleeding surface using some pressure. The period of time pressure is needed and the amount of INSTAT MCH necessary to achieve hemostasis will be dependent upon the nature and amount of bleeding to be controlled. It has been shown that hemostasis usually occurs within 2-4 minutes. In contact with blood, INSTAT MCH becomes a coherent gelatinous unit that is completely absorbed if left *in situ*. However, it is suggested that the surgeon simply remove the collagen mass prior to wound closure.

#### **ADVERSE REACTIONS**

Adverse reactions reported with another microfibrillar collagen hemostatic agent that were possibly related to its use were adhesion formation, allergic reaction, foreign body reaction, and subgaleal seroma (report of a single case). The use of microfibrillar collagen in dental extraction sockets has been reported to increase the incidence of alveolalgia.

#### **HOW SUPPLIED**

Sterile INSTAT MCH, an absorbable collagen hemostatic agent in a microfibrillar form, is supplied in the following sizes:

Product Code Size Quantity Per Case

1985 1 gram 12

1984 0.5 gram 12

Do not resterilize this product.

#### **CAUTION**

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

## ***Thrombin, Topical USP THROMBOGEN (Bovine Origin)***

### **THROMBIN, TOPICAL USP**

### **THROMBOGEN (BOVINE ORIGIN)**

THROMBOGEN Thrombin must not be injected! Apply on the surface of bleeding tissue as a solution or powder.

### **DESCRIPTION**

(THROMBIN, TOPICAL USP) THROMBOGEN is a protein substance produced through the conversion reaction in which prothrombin of bovine origin is activated by tissue thromboplastin in the presence of calcium chloride. It is supplied as a sterile powder that has been freeze-dried in the final container. Also in this preparation are calcium chloride, sodium chloride, aminoacetic acid (glycine) and benzethonium chloride. Glycine is included to make the dried product friable and more readily soluble, and benzethonium chloride as a preservative at 0.2 mg./vial. This product is prepared under rigid assay control. The unit is defined by the U.S. standard and is approximately equal to the stated unitage. A U.S. unit is defined as the amount required to clot 1 ml of standardized fibrinogen solution in 15 seconds. Approximately 2 U.S. units are required to clot 1 ml of oxalated human plasma in the same period of time.

### **CLINICAL PHARMACOLOGY**

THROMBOGEN Thrombin requires no intermediate physiological agent for its action. It clots the fibrinogen of the blood directly. Failure to clot blood occurs in the rare cases where the primary clotting defect is the absence of fibrinogen itself. The speed with which thrombin clots blood is dependent upon its concentration.

### **INDICATIONS AND USAGE**

THROMBOGEN Thrombin is indicated as an aid in hemostasis wherever oozing blood from capillaries and small venules is accessible.

In various types of surgery, solutions of THROMBOGEN Thrombin may be used in conjunction with Absorbable Gelatin Sponge, USP for hemostasis.

### **CONTRAINDICATIONS**

THROMBOGEN Thrombin is contraindicated in persons known to be sensitive to any of its components and/or to material of bovine origin.

### **WARNING**

Because of its action in the clotting mechanism, THROMBOGEN Thrombin must not be injected or otherwise allowed to enter large blood vessels. Extensive intravascular clotting and even death may result. THROMBOGEN Thrombin is an antigenic substance and has caused sensitivity and allergic reactions when injected into animals.

### **PRECAUTIONS**

General--Consult the Absorbable Gelatin Sponge, USP product labeling for complete information for use prior to utilizing the thrombin-saturated sponge procedure.

Pregnancy--Category C--Animal reproduction studies have not been conducted with Thrombin, Topical USP (Bovine origin). It is also not known whether Thrombin, Topical USP (Bovine origin) can cause fetal harm when administered to a pregnant woman or can affect reproduction capacity. Thrombin, Topical USP (Bovine origin) should be given to a pregnant woman only if clearly indicated.

Pediatric Use--Safety and effectiveness in children have not been established.

### **ADVERSE REACTIONS**

Allergic reactions may be encountered in persons known to be sensitive to bovine materials.

### **DOSAGE AND ADMINISTRATION**

Solutions of THROMBOGEN Thrombin may be prepared in sterile distilled water or isotonic saline. The intended use determines the strength of the solution to prepare. For general use in plastic surgery, dental extractions, skin grafting, neurosurgery, etc. solutions containing approximately 100 units per ml are frequently used. For this, 10 ml of diluent added to the 1,000 unit package is suitable. Where bleeding is profuse, as from abraded surfaces of liver and spleen, concentrations as high as 1,000 to 2,000 units per ml may be required. Intermediate strengths to suit the needs of the case may be prepared by selecting the proper strength package and dissolving the contents in an appropriate volume of diluent. In many situations, it may be advantageous to use THROMBOGEN Thrombin in dry form on oozing surfaces.

**Caution:** Solutions should be used immediately upon reconstitution. However, the solution may be refrigerated at 2-8° C for up to three hours.

The following techniques are suggested for the topical application of THROMBOGEN Thrombin.

1. The recipient surface should be sponged (not wiped) free of blood BEFORE THROMBOGEN THROMBIN IS APPLIED.
2. A spray may be used or the surface may be flooded using a sterile syringe and small gauge needle. The most effective hemostasis results when the THROMBOGEN Thrombin mixes freely with the blood as soon as it reaches the surface.
3. In instances where a concentration of approximately 1,000 units per ml is desired, the contents of the vial of sterile isotonic saline diluent may be easily transferred into the vial of THROMBOGEN Thrombin with a sterile syringe or a sterile transfer needle. If the transfer needle method is utilized, transfer the diluent as follows:
  - (a) Flip plastic cover off of diluent vial.
  - (b) Remove clear plastic cover from the transfer needle by twisting to break the seal. (c) Insert the exposed needle into the diluent.
  - (d) Flip the plastic cover up on the THROMBOGEN Thrombin vial. **DO NOT REMOVE THE COVER AND ALUMINUM SEAL.**
  - (e) Remove the pink plastic cap from the transfer needle exposing the needle.
  - (f) Invert the vial of diluent and insert the exposed needle into the THROMBOGEN Thrombin vial.
  - (g) Allow the vacuum to draw the complete contents of the diluent vial into the THROMBOGEN Thrombin vial.

#### **IMPORTANT**

\* Insert needle in diluent first.

\* Do not remove the cover and aluminum seal on the THROMBOGEN Thrombin vial

4. In instances where THROMBOGEN Thrombin in dry form is needed, the vial is opened by removing the metal ring by flipping up the metal tab and tearing counterclockwise. The rubber-diaphragm cap may be easily removed and the dried THROMBOGEN Thrombin is then broken up into a powder by means of sterile glass rod or other suitable sterile instrument.

5. Sponging of treated surfaces should be avoided in order that the clot remains securely in place.

THROMBOGEN Thrombin may be used in conjunction with Absorbable Gelatin Sponge, USP as follows:

1. Prepare THROMBOGEN Thrombin solution of the desired strength.
2. Immerse sponge strips of the desired size in the THROMBOGEN Thrombin solution. Knead the sponge strips vigorously with moistened fingers to remove trapped air, thereby facilitating saturation of the sponge.
3. Apply saturated sponge to the bleeding area. Hold in place for 10 to 15 seconds with a pledget of cotton or a small gauze pad.

#### **HOW SUPPLIED**

THROMBOGEN Thrombin is supplied as:

NDC 56091-038-01 Package contains one 1,000 unit vial of THROMBIN, TOPICAL USP THROMBOGEN.

NDC 56091-039-05 Package contains one 5,000 unit vial of THROMBIN, TOPICAL USP THROMBOGEN with one 5 ml vial of isotonic saline diluent and a transfer needle.

NDC 56091-040-05 Package contains one 10,000 unit vial of THROMBIN, TOPICAL USP THROMBOGEN with one 10 ml vial of isotonic saline diluent and a transfer needle.

NDC 56091-041-05 Package contains one 20,000 unit vial of THROMBIN, TOPICAL USP THROMBOGEN with one 20 ml vial of isotonic saline diluent and a transfer needle.

#### **STORAGE**

Store at normal, controlled room temperature, except 1,000 unit vial of THROMBIN, TOPICAL USP THROMBOGEN must be stored between 2?-8?C (36?-46?F).

## ***POWERSTAR Bipolar Scissors***

### **POWERSTAR BIPOLAR CABLE**

#### **INDICATIONS**

The POWERSTAR\* Bipolar Scissors are non-sterile, reusable devices intended to facilitate cutting and bipolar coagulation of soft tissue in open surgical procedures.

#### **CONTRAINDICATIONS**

The POWERSTAR Bipolar Scissors are not indicated for contraceptive coagulation of the fallopian tube, but may be used to achieve hemostasis following transection of the fallopian tube.

#### **WARNING AND CAUTIONS**

Do not use this device with a monopolar source output; damage to the instrument and harm to the patient or medical personnel could result.

A thorough understanding of the principles and techniques associated with electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical isolation or grounding is not compromised.

Avoid contact of the bipolar scissors with a monopolar instrument. Damage to the electrosurgical generator may result.

Limit the use of the scissors to only properly trained, qualified medical personnel.

Because of the variability of output voltages and modes from generator to generator, do not use the instrument with generators having bipolar output voltages that exceed 1000 volts peak to peak. Refer to appropriate electrosurgical generator manual for indications and instructions on bipolar output characteristics to ensure that all safety precautions are followed.

The POWERSTAR Bipolar Scissors are only compatible with the POWERSTAR Bipolar Cable. Use of these scissors with other cables will result in damage to the scissors.

To initiate the use of the POWERSTAR Bipolar Scissors, the power setting on the generator should be set to the lowest power setting needed to achieve adequate hemostasis. If adequate hemostasis is not achieved, increase power settings in small increments until adequate hemostasis is achieved.

Never touch the active blades of the scissors while energized. This may result in burns. Use surgical gloves designed by their manufacturers for electrosurgical procedures. The use of other gloves may result in burns or shocks.

Do not touch the blades of the Bipolar Scissors to any staples, clips, or sutures when the scissors are energized. Damage to the Bipolar Scissors and these products may occur.

Do not use the POWERSTAR Bipolar Scissors to cut staples or clips. When not energized, the scissors can be used in a manner appropriate for standard surgical dissecting scissors of like type and size. Cutting of sutures, meshes and similar materials may lead to premature dulling of the scissors blades.

Scissors are not intended to be sharpened or repaired in any way. This may damage the electrodes.

Do not use abrasive cleaning tools or solutions on the device.

#### **CAUTION**

Federal law (U.S.) restricts this device to sale by or on the order of a physician.